

AE-456: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h

Date: 2025-03-04


Tags: Dihydroxo Ag₂O PNP
P(tBu)N(py)P(tBu) $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$
NMR ATR-IR AE Future reference
analytics ¹H ³¹P IR Elemental Analysis

Category: Two-photon

Status: Done

Created by: Alexander Eith

Reaction scheme/sample structure



Literature/reference experiments

Literature	https://doi.org/10.1039/D1EE01053K
Reproduction	/
Related experiments	Experiment - AE-KRA-204: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h Experiment - AE-238: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h Experiment - AE-252: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h AE-310: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h Two-photon - AE-382: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h

Reagents

Name	CAS-Number/Experiment-Number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume _{theo} [ml]	Volume _{exp} [ml]	Density [g/mL]
$[\text{RuCl}_2(\text{CO})(\text{PNP}\{\text{tBu}\})]$	Two-photon - KRA-DS-02: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$ I	0.42	1.00	250.00	248.84	595.53	/	/	/
Ag_2O	20667-12-3	2.58	6.14	597.32	596.73	231.74	/	/	/
H_2O	Prep work - AE-442: Degassing of Milli-Q water	98.23	234	1770.13	/	18.02	1.774	1.8	0.998
THF	Prep work - AE-396: Degassing of THF	/	/	/	/	/	12.5	13	/

Work-up Reagents

Name	CAS Number / Experiment-Number	Inventory number	Mass _{exp} [g]	Volume [ml]	Concentration [M]
THF	Prep work - AE-396: Degassing of THF	/	/	13.5	/
n-heptane	AE-315: Degassing of n-Heptane	/	/	20 + 4.5	/

Procedure/observations

All steps, unless mentioned otherwise, were carried out under argon using standard Schlenk technique or in an argon filled glovebox. (see Protocol [Protocol - Schlenk Technique](#))

Due to the instability of Ag₂O and Dihydroxo to the light, the procedure was done mostly in the dark, by either wrapping the flasks containing either or both in Aluminium foil or by tuning off the light in the lab and fumehood

Date	Time	Step	Photos	Observations
04.03		A 25 mL brown glas Schlenk flask was prepared, by heating under vacuum		
	9:35	[RuCl ₂ (CO)(PNP)] and Ag ₂ O was weighted in air into the reaction schlenk flask and the flask was set under argon.		
		The light in the fumehood was turned out		
	10:04	THF (13 mL) was added to the reaction schlenk flask.		
	10:05	Water (1.8 mL) was added to the reaction schlenk flask		
	10:05 - 15:05	The reaction mixture was stirred for 5 h at 550 rpm. During this the flask was wrapped with aluminium foil		
	15:05 - 15:07	The reaction mixture was allowed to settle		
	15:10 - 15:15	The mixture was filtered through a pore 3 frit according to Protocol - Filtration with frit technique into a 100 mL-Schlenk	1st dilter residue.jpg 1st filtrate.jpg	Significant amounts of the solid went through the filter frit, dark solid, purple/black filterate
		A new 100 mL flask was prepared		
	15:35 - 15:42	The obtained suspension was filtered through a PTFE syringe filter (pore diameter: 0.2 µM) in the new 100 mL flask	after filtration through PTFE filter.jpg	dark violet, clear solution

	15:47 - 15:55	The solvents of the filtrate were removed under reduced pressure.		approx. 2.5 mL remained in the flask
		The water bath was removed		
	15:55 - 16:08	The remaining solvents were removed without water bath, trying to freeze the solution.	after removing water bath – bad quality.jpg	Did not work, remained liquid. Bubbeled up during drying, still the solvents were removed and a solid was obtained Purple suspension
	16:08 - 16:15	The obtained solid was dried under reduced pressure	after removing solvents.jpg	Pruble brownish film
		The solid was stored in the dark at rt		
05.03	11:00 - 11:15	The obtained solid was dried under reduced pressure		
	11:17	To the solid THF (13.5 mL) was added und the obtained mixture was shaken strongly for approx. 2 min and than stirred for approx. 8 min at 665 rpm	in THF.jpg	purple brwonish suspension
	11:27	n-heptane (20 mL) was added to the mixture	after addition of heptane.jpg	purple solution, brownish solid
	11:30 - 15:00	The mixture was filtered through a pore 4 frit according to Protocol - Filtration with frit technique into a 100 mL-Schlenk		very slow filtration, but still working.
	15:00 - 15:07	The solid was washed with n-heptane (3 * 1.5 mL)	2nd filtrate.jpg	purple solution
	15:10 - 17:30	The solid was dried under reduced pressure	product.jpg	
		The solid was stored in Equipment - Freezer Lab E004 CEEC II at -20 °C in the dark under Argon		
11.03	9:00	The solid was transffered into a 2 mL vial in a 10 mL Schlenk flask		Some material was lossed, since it sticked stongly to any surface
		During this an aliqoute (approx. 2 mg) was transferred into another 10 mL Schlenk flask		
	9:15	Both flask were set under argon		

	13:35	From the aliquote n NMR sample was prepared in MeCN-d3 (0.7 mL) according to Protocol - Preparation of NMR Sample	NMR sample.jpg	AE-456-1 yellowish solution
01.04		5 mL NS 14 schlenk flask was prepared		
	10:10	To the flask AE-456-1 (approx. 15 mg) was added		
		The flask was set under argon		
	10:33	To the flask MeCN-d3 (0.7 mL) was added		
	10:34	The flask was stirred at 500 rpm	dissolved.jpg	
	10:38	The NMR was prepared according to Protocol - Preparation of NMR Sample	in NMR.jpg	AE-456-1-1
02.04	15:20	AE-456-1-1 was submitted	before NMR.jpg	
04.06	09:40	IR spectra where mesured according to Protocol - IR measurements, Lab E007, CEEC I using the Equipment - Fourier Transform Infrared Spectrophotometer Shimadzu IRSpirit Several conditions for the measurment where tried: * AE-456-1-IR-bkg: background measurment * AE-456-1-IR-1: resolution: 4 cm ⁻¹ , Absorbtion mode, 30 scans * AE-456-1-IR-2: resolution: 4 cm ⁻¹ , Absorbtion mode, 30 scans * AE-456-1-IR-3: resolution: 8 cm ⁻¹ , Absorbtion mode, 30 scans * AE-456-1-IR-4: resolution: 8 cm ⁻¹ , Absorbtion mode, 30 scans * AE-456-1-IR-5: resolution: 4 cm ⁻¹ , Absorbtion mode, 30 scans * AE-456-1-IR-6: resolution: 16 cm ⁻¹ , Absorbtion mode, 30 scans		
21.07	11:20	Two flask (flask 1 and 2) were prepared equipped each with a 2 mL screw cap vial		
	11:40	To flask 1 AE-456-1 (3.5 mg) was added		AE-456-EA
	11:45	To flask 2 AE-456-1 (approx. 4 mg) was added		AE-456-IR-a
	13:20	From AE-456-IR-a IR spectrums were meausred in air by Akuila Edwards (AG Schmittty)		
	14:00	AE-456-EA was submitted for elemental analysis		
22.07	/	The elemental analysis of AE-456-EA was measured by Sandra Köhn		

Analysis

Date	Time	Sample name	Analysis method	Used device	Solvent	Raw Data	Processed Data	Comparative Data	Interpretation
12.03	11:12	AE-456-1	NMR 1H, 31P{1H}quant	IAAC 400 MHz I	MeCN-d3	AE-456-1.zip	AE-456-1_10.nmrium	Two-photon - AE-382: Synthesis of [Ru(CO)(OH)2(PNPtBu)] using Ag2O, 5 h (AE-382-1)	As expected. Purer than AE-382, very few impurities, also in 1H

02.04	15:26	AE-456-1-1	NMR 1H, 13C{1H}, 31P{1H}, 13C{1H}-APT, 1H/13C-COSY, 1H/13C{1H}-HMBC, 1H/13C{1H}-HSQC, 1H/31P{1H}-HMBC	1000 500 MHz	MeCN-d3	AE-456-1-1-1H.zip AE-456-1-1-13C.zip AE-456-1-1-31P.zip AE-456-1-1-31P-HMBC.zip AE-456-1-1-HMBC.zip AE-456-1-1-HSQC.zip AE-456-1-1-COSY.zip AE-456-1-1-13C-APT.zip	AE-456-1-1-1H-13C-31P-APT.nmrium AE-456-1-1-31P-HMBC.nmrium AE-456-1-1-COSY.nmrium AE-456-1-1-HMBC-HSQC.nmrium	AE-456-1	Full evaluation pending. 1H and 31P fit with previous 13C: 2 weak triplets for CH3 groups, 3 triplets for Cq(CH3) and CH2, 4 signals in aromatic area, 3 expected, 1 signal for CO (very weak) 31P HMBC: as expected, no other signals observed 13C HSQC: CH3 and CH2 groups observed, aromatic as well 13C HMBC: CH3/CH2 coupling observed, CH2/Cq(CH3) ?, CH2/Cq(Ar) observed, no other coupling apparent in Aromatic are at first glance COSY: Aromatics with each other, CH2 with each other, CH3 with itself
04.06	09:40	AE-456-1-IR	ATR-IR	Protocol - IR measurements, Lab E007, CEEC I	/	AE-456-1-IR-bkg.txt AE-456-1-IR-1.txt AE-456-1-IR-2.txt AE-456-1-IR-3.txt AE-456-1-IR-4.txt AE-456-1-IR-6.txt AE-456-1-IR-2.ispd AE-456-1-IR-3.ispd AE-456-1-IR-4.ispd AE-456-1-IR-5.ispd AE-456-1-IR-6.ispd	Used for plotting: plot_ir_spectra.py AE-456-1-IR-all.png AE-456-1-IR-transmission.png AE-456-1-IR-good_spectra.png Used for plotting and evaluation: IR_single_spectra.py AE-456-1-IR-3.png AE-456-1-IR-3_peaks.csv AE-456-1-IR-4.png AE-456-1-IR-good_spectra.png	/	In absorbance mode no good spectra was obtained In transmission mode good spectra were obtained only with a resolution of 8 cm ⁻¹ , with higher and lower resolution artefacts (like positive peaks where observed) -3 and -4 are very similar (also same experimental conditions). For them full evaluation was done. CO band can be clearly seen, in water area some distinct peaks. In general ok signal strength (down to 65 % transmission).
21.07	13:20	AE-456-IR-a (Name in file: Eith_a)	ATR-IR with diamond	IR from AG Schmitt (measured by akulla)	/	Eith_a.CSV	Used for plotting: IR_new.py Eith_a_peaks.csv AE-456-IR_Eith_a.png	AE-456-1-IR Used background: Background-Mon-Jul-21-13-17-01-2025-GMT02-00.CSV	Better signal to noise ratio, still quite weak signals In general very similar to AE-456-1-IR
21.07	13:20	AE-456-IR-a (Name in file: Eith_a1_200scans)	ATR-IR with diamond	IR from AG Schmitt (measured by akulla)	/	Eith_a1_200scans.CSV	Used for plotting: IR_new.py AE-456-IR_Eith_a1_200scans.png Eith_a1_200scans_peaks.csv	AE-456-1-IR Used background: Background-Mon-Jul-21-13-17-01-2025-GMT02-00.CSV	Better signal to noise ratio, still quite weak signals Best quality IR so far In general very similar to AE-456-1-IR
22.07	/	AE-456-EA (Name in file: AE-456EA)	elemental analysis	elemental analysis from AG Schubert (measured by Sandra Köhn)	/	AE-456-EA.jpg	C: 50.48 % H: 7.13 % N: 2.68 %	C: 51.6 % H: 8.12 % N: 2.51 % With additional 1 H2O: C: 49.99 % H: 8.22 % N: 2.43 % With additional CO2 - H2O C: 51.36 % H: 7.41 % N: 2.40 % Further possible molecules: AE-456-EA-calc.xlsx Old data: Two-photon - AE-382: Synthesis of [Ru(CO)(OH)2(PNPtBu)] using Ag2O, 5 h	Significant error compared to theoretical composition, especially H has to large deviation. Can also not be explained by addition of not carbon containing molecules to lower H content in theoretical values. But at least significant better than AE-382-EA Other impurities cannot help, free ligand and coordinated/crystallized MeCN, THF, n-heptane all would yield increased C and H content. Inorganic compounds (like Ag2O) do help, with 2.5 % we get C: 50.31 %, H: 7.92 % and N: 2.45 %, C and N good, but H still to high, more Ag2O: unrealistic and yields to low C and N No realistic impurity found to explain deviations in elemental analysis. Maybe error in measurement, but unlikely Errors up to 1 % can realistic appear for pure, simple organic compounds: https://doi.org/10.1021/acscentsci.2c00325 , just by random error, will be larger if sample is not 100 % pure and more deficit like here

Product characterization and results

Sample	Mass [mg]	Purity	Mass _{pure} [mg]	Amount [μmol]	Yield [%]	Description	Storage location
AE-456-1	91.32	> 97 %	/	0.163	39	greenish grayish solid, product.jpg	Equipment - Freezer Lab E004 CEEC II A

Future recommendations

Old procedure	Problem	Suggested new procedure
1st filt through pore 3	not all solids removed	use pore 4 or 5

Linked experiments

- [AEI-010: Synthesis of \[Ru\(CO\)\(OH\)2\(PNN\)\] using Ag2O](#)

- JSC-KS-04: Degassing of n-Heptane
- AEI-024: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{NNN})]$ with 2,6-bis(N,N-diethylaminomethyl) using Ag_2O
- AEI-026: Capillary production
- AEI-048: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{NNN})]$ with 2,6-bis(N,N-diisopropylaminomethyl) using Ag_2O
- AEI-064: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$ with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag_2O
- AEI-073: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$ with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag_2O
- AE-101: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$
- AE-102: Degassing of THF
- AE-105: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$ with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag_2O
- AE-108: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$
- AE-114: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$ with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag_2O
- AE-121: Degassing of n-Heptane
- AE-145: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$ with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag_2O
- AE-152: Degassing of MeCN
- AE-157: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$
- AE-167: Irradiation of PhPDA (365 nm + white LED, gas phase measurement), AE-151, 1.5 mg/mg SDS, 2 mg/mL PhPDA
- AE-198: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP}\{\text{tBu}\})]$ with Ag_2O with different reaction times (1, 2, 5, 24 h)
- AE-KRA-204: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h
- AE-208: Degassing of Milli-Q water
- KRA-080: Degassing of THF
- KRA-104: Degassing of Milli-Q water
- AE-230: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via

[RuCl₂(CO)(p-cymene)]

- AE-232: Synthesis of [RuCl₂(CO)(PNP)] with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via [RuCl₂(CO)(p-cymene)]

- AE-238: Synthesis of [Ru(CO)(OH)₂(PNPtBu)] using Ag₂O, 5 h

- AE-246: Radiation of [Ru(CO)(OH)₂(PNP{tBu})] (AE-238, 97 % pure) 365 nm and white

- AE-252: Synthesis of [Ru(CO)(OH)₂(PNPtBu)] using Ag₂O, 5 h

- AE-261: Synthesis of [RuCl₂(CO)(PNP)] with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via [RuCl₂(CO)(p-cymene)]

- AE-281: Stability test of [Ru(CO)(OH)₂(PNP{tBu})] (AE-252) under air and afterwards irradiation with 365 nm and white

- AE-310: Synthesis of [Ru(CO)(OH)₂(PNPtBu)] using Ag₂O, 5 h

- AE-315: Degassing of n-Heptane

Prep work - AE-396: Degassing of THF

Prep work - AE-442: Degassing of Milli-Q water

Two-photon - AE-382: Synthesis of [Ru(CO)(OH)₂(PNPtBu)] using Ag₂O, 5 h

Two-photon - KRA-DS-02: Synthesis of [RuCl₂(CO)(PNP)] with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via [RuCl₂(CO)(p-cymene)] I

Linked resources

Equipment - [Fourier Transform Infrared Spectrophotometer Shimadzu IRSpirit](#)

Equipment - [Freezer Lab E004 CEEC II](#)

Protocol - [Preparation of NMR Sample](#)

Protocol - [Filtration with frit technique](#)

Protocol - [Schlenk Technique](#)

Protocol - [Preparation of capillary for H₂O-NMR](#)

Protocol - [IR measurements, Lab E007, CEEC I](#)

Attached files

AE-456-EA-calc.xlsx

sha256: 9c08760ff858b5fb761784f31a62ebfddc22cbce4fda2f5afc2024076d921a81

AE-456-EA.jpg

sha256: aceb21b2125bf4562df176b9067bf6ed868267a863629b6d720811b7022da41d




Background-Mon-Jul-21-13-17-01-2025-GMT02-00.CSV

sha256: 9039f96dff4009c94f187d0cf2637100d2d437187b4f97fb33b7e7ea0d1b7cb8

AE-456-IR_Eith_a.png

sha256: 5e8c2a78bdd8737b9b6928a48df00092bbc23c4a60dea88e8c34788e15373f1c



IR_new.py


sha256: 5f67422508287e2ed7aeacd4c38e2cd4f8948737194f4fc45fa0206c90f2a086

Eith_a_peaks.csv

sha256: e27415521dd56c8ffe4113c9638b337f1f3892e3ed4fb130839c30e2af16bb4c

AE-456-IR_Eith_a1_200scans.png

sha256: 55ff7f233a4cda7308c8264058445014d6f78e3f6a158fa9dcd9a2eddaae03fe



Eith_a1_200scans_peaks.csv

sha256: abd9627ac1e1de04776431e38dee3cdae59c2af7ff70f350c644d0b7f2b0eb90

Eith_a.CSV

sha256: 3907b1a1e8c19ddd8413f8008cd4b6575152f754d1dd390b6c36bb7674dfa24b

Eith_a1_200scans.CSV

sha256: 3f0721ed61420fe21f227c8c275adbd6fbe324ed00cac61a98a8a06dd390ff6d

AE-456-UVVis_2_7420287SP.TXT

sha256: 0ac957f02f10e345c9227a3126fa521d463dd92f9234177403705461354ad85b

AE-456-UVVis_2.ABS8

sha256: 5b9317bee247a1cd3fb4d80cc750c29bfc70edc3fafbb2934f2fed5e672812e

AE-456-UVVis_1_7420287SP.TXT


sha256: bc6c8be7c224d1f382b0c82c819c10ca0d67658afa2bdba5345d0a330bb741e3

AE-456-UVVis_1.ABS8

sha256: f8b0baadd3c0c73b835d4fda7fcf7744da1665a2ce3c006feb97e2c0266c2e91


AE-456-1-IR-all.png

sha256: 18e0c999ba9de1c4ee52fc581b08cd63ea5b9fc55a07a3319d2f1cc8d1dcbf62




AE-456-1-IR-4.png

sha256: 32442f8b3814e1e4e11034e8c479f10f86b451aa3dcfe98702c1ef6b58b1e7f9




AE-456-1-IR-3.png

sha256: ceae002988bbbe66e54796bb0b2b9143ee7aa8bbeaa18254c526756e3ae0aadd



AE-456-1-IR-good_spectra.png

sha256: fa6bf361c3d705bb6fc565001b94d8ed71c981a7356e8c6059a21b5c9250e728



AE-456-1-IR-transmission.png

sha256: e88c4912e56b8aa59226c0e62795821ef0da61a69101dff25a12877244b007b4



AE-456-1-IR-4_peaks.csv

sha256: 2846f57eb43dfafed6f4dfe02c68f12c24d730ff261ffc62d245f7fe96254a57

AE-456-1-IR-3_peaks.csv

sha256: 9eb9b49cd7466f3ad8f54dab798535eb710d57b467968db76979211c204b32cb

IR_single_spectra.py

sha256: 6bc4370ea64168ed4fc79c0ca5fa65bebb76bdc8b587309344c00693e6c58c97

plot_ir_spectra.py

sha256: fd415843ff163d159d56ac3f3c160db712cb7d140f4dceee7874f8732fe22d2e

AE-456-1-IR-2.txt

sha256: dbdefd894ac58ab395d8a32a600e8047137e2db60ad760b40e2d17e84cc00205

AE-456-1-IR-6.txt

sha256: 85f042f5d9607064864c3cbb995df42e5482c1dd11f3b1aa469e5938b07a8e29

AE-456-1-IR-6.ispd

sha256: c827f18f8de5204b6776d2a58aae2496e4faf83edfbeb6f8334f9a1196b43700

AE-456-1-IR-5.ispd

sha256: 852350eb9bccead0ffa7557f9cf6b84e1e3ade2f453f09d569aa79b762fc3202

AE-456-1-IR-4.txt

sha256: c5b568061fedd33abb3fb57811ce1683fd96b9c12a41f03900a0b3b7e77b25bf

AE-456-1-IR-4.ispd

sha256: 5aa01885aa5378601dcbf87e0e3481cfdbdb55a1606dfcaa2da9522d347d18af

AE-456-1-IR-3.txt

sha256: 6233dd6d6d95ab913d9d5123e7a3866ad0585150b4c4bee0d87ebdcfc6bec61a

AE-456-1-IR-3.ispd

sha256: 0a2a89827e50cae3eba94d0656a070669063bd843a445ce36f1ad0b311726df8

AE-456-1-IR-2.ispd

sha256: 4d6ed27c532de15109128e8551efb1f6ec595c611a1f4d67cdeb5fd262d0a06d

AE-456-1-IR-bkg.txt

sha256: 5d854d047fe9fc8a21a88b82743a4fb53115e3a76160541d4bb523e774425337

AE-456-1-IR-1.txt

sha256: 93d4129204d7fe3ae75a41f86b1515dc5626b02854205198f40a68d57ad22491

AE-456-1-1-HMBC-HSQC.nmrium

sha256: bd644af3cccb0e2345e7c1d9e8f72aad9a3b61dbfd4562039402a6e0e16d3979

AE-456-1-1-COSY.nmrium

sha256: 64eb40edb89b79397d78764eb7d2924e5cd5e45d4623f0305b9c79ef34da171e

AE-456-1-1-31P-HMBC.nmrium

sha256: 4c231487eedea1cf9dd1e2c436c9cc466817a32d65654ebb739e09172b21cb8e

AE-456-1-1-1H-13C-31P-APT.nmrium

sha256: 4c6613cd0692ecd7411931f28a596ef5746587cd5f7b794dd08c2b64adc09c7a

AE-456-1-1-13C-APT.zip

sha256: b75622a2c584c9c7d51baf9c5c221babe7b9f2264d9701398a7f3cac570af50c

AE-456-1-1-COSY.zip

sha256: 652358325b7c3f60ab13a7d64c39c2ee98b66e36b1a4379572476c20bf3fbd47

AE-456-1-1-HSQC.zip

sha256: 559c1b413bb4bb6e355f04811c462974b49205d2fb580a9fd011c30859e697dc

AE-456-1-1-HMBC.zip

sha256: b387c112008a3293b1ec162f6eddb5f7479a93c7d43455eede667ef1c4a72f3e

AE-456-1-1-31P-HMBC.zip

sha256: a8c65071ac5964a5ff645c227a4fbd68b83632a816d1b4dcba23eb247a817ed9

AE-456-1-1-31P.zip

sha256: 51579b7704070c3e70a0699a8f86b510d424dd2620c6cea0afa5c4bffec4ed13

AE-456-1-1-13C.zip


sha256: bff1c5e084731179249ab8b7af203bbb59467cc30bcbd87badcd3effcf79e4a4

AE-456-1-1-1H.zip

sha256: 571fd463f56687e21bda334f3ab90c4da037c89ad0628640409cc2a97f1074c3

before-NMR.jpg

sha256: 189d713f520f318e615ce478c412b71fca49dfcf2013fd4765f138b947d5e19e



dissolved.jpg

sha256: 8878cac4a231e94d984ae3117353b15155ca0d1377550f23674ede27ecbc18a9




in-NMR.jpg

sha256: 806e30c941773c544fe098a300da23991d21ec7403efbd9d86e3fe7e2232cb36




after-filtration-through-PTFE-filter.jpg

sha256: 37e8cf496a23e55b57a9006196b8cdcb987034a3839b9372a6eb0c05ba097bc




NMR-sample.jpg

sha256: bda6452ea894e0dd30d7dcd6b16d68dbf265c964124162ba3379c149cd29df4b



2nd-filtrate.jpg

sha256: c07a8ac85aeeb76636e78c80862299858239caeeb5c4a1c1abae2fe803d53d1f



product.jpg

sha256: e15f16e125f8198f33098a7c968e96c52eda190f247f30e2fdf3b77d6c268d88




after-addition-of-heptane.jpg

sha256: 42c88ac228a5298b5479164af0845240710b9bed34e9512f1020986462ff7e6e



after-removing-solvents.jpg

sha256: 17fe480f13e4a6953d8334fcc0038e42b47fc4889c55a7ca6da894928fdc817e



in-THF.jpg

sha256: 4feabb176157302c940597b88ed802b8d0369efcb1412d941a558ef9fe2e37ae



1st-filtrate.jpg

sha256: 7283aab0a48c4890cd84cb733fc8e2be7f747a262effd880e0bf79448416bcc1




1st-dilter-residue.jpg

sha256: 1c411f28cb1c470e3debe8e9645ff5cfe87b1d2190e446b335f45c88bc5a0905



after-removing-water-bath-_-bad-quality.jpg

sha256: 85a3d91b7e5f8afde9acfe88077c31df5e6c677f5a4224b5247422eba4ed824d



AE-456-1.zip

sha256: 464fa58808dfdeeff5c479df021a1e2a4f4db5fab49d58b2744448ae630ab467

AE-456-1_10.nmrium

sha256: 2f90c518948b7e8712bf976a9e2c0d4ec8b2be257b37246f1070b807c74b1ee8

AE-456.cdxml

sha256: 70d2a67f6ae7477a3c31ec8e776e001181c617b18e0d276cd2995bdbdfdb2ca2



Unique eLabID: 20250304-419a10d9bb233002c38616635e41b7417c9d8e27

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1825>