

AE-441: Determination of absorption coefficient of [Ru(CO)(OH)2(PNP)] with UV/Vis

Date: 2025-02-05
Tags: Dihydroxo PNP P(tBu)N(py)P(tBu) [Ru(CO)(OH)2(PNP)] AE UV/Vis
Category: Two-photon
Status: Done
Created by: Alexander Eith

Reaction scheme/sample structure

Literature/reference experiments

Literature	https://doi.org/10.1039/D1EE01053K
Reproduction	/
Similar experiments	PBDF - NB-JSC-143: Testing of BFD2 (MZ-124-6) DMSO and DMSO-H2O solutions

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [μmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]
[Ru(CO)(OH)2(PNP)]	Two-photon - AE-382: Synthesis of [Ru(CO)(OH)2(PNPtBu)] using Ag2O, 5 h (AE-382-3)	/	4.47 + 4.0	/	2.5 + 2.233	2.58 + 2.22	558.65	/	/
water	Prep work - AE-442: Degassing of Milli-Q water	/	/	/	/	/	/	/	2.5 + 4.0 + 1.5 + 1.4 + 1.5 + 1.8

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard [Protocol - Schlenk Technique](#) or in an nitrogen filled glovebox (MBRAUN).

The synthesis/reactions were performed in the dark by either wrapping the reaction flasks in aluminum foil or turning the lab light and light in the fume hood off, when the flasks were uncovered.

Date	Time	Step	Pictures	Observations
05.02		2 10 mL Schlenk flask were set under argon		
	10:10	To flask 1 [Ru] (2.58 mg) and to flask 2 [Ru] (2.22 mg) were added in air		
	10:20	Both flasks were set under argon: flask 1: AE-441-1 , flask 2: AE-441-2		

	11:53 - 11:56	To -1 water (2.5 mL) and to -2 water (4.0 mL) was added		
	- 12:16	The obtained solutions were stirred at 500 rpm	left: -2, right -1.jpg	left: clear, yellow solution, right: cloudy, yellow-orange suspension
		The Equipment - Avantes UV/vis spectrophotometer was turned on and allow to warm up for approx. 10 min		
		Both lamps were tuned on and the device was connected to the computer		
		A dark refernce spectra was measured to assure, that the light intesitiy is fitting		
	12:12	A reference spectra of water was measured.	water.jpg	clear
		The water was discarded and the cuvette was dried by shaking off most of the water		
		From here on AE-441-1		
		The suspension was taken into a 5 mL syringe		
		From now work under air for -1:		
		The suspension was filtered through a PA syringe filter (pore diameter = 0.2 µm) into the quartz cuvette used for the reference measurement of water.		
	12:21	An UV/Vis was recorded AE-441-1-1	-1-1.jpg	To concentrated
		The sample was pured in a 50 mL beaker		
		The cuvette was washed once with water and dried by shaking off most of the water		
		from AE-441-1-1 0.5 mL were taken off and transferrred into the cuvette using a 1 mL Eppendorf syringe		
		To the cuvette water (1.5 mL) was added using a 1 mL Eppendorf syringe		
	12:25	An UV/Vis was recorded AE-441-1-2	-1-2.jpg	Better, still to concentrated
		The sample was pured in a 10 mL snap on cap vial		
		The cuvette was washed once with water and dried by shaking off most of the water		

		from AE-441-1-2 0.7 mL were taken off and transferred into the cuvette using a 1 mL Eppendorf syringe		
		To the cuvette water (1.4 mL) was added using a 1 mL Eppendorf syringe		
	12:28	An UV/Vis was recorded AE-441-1-3		Good intensity
		The sample was pured in a 10 mL snap on cap vial		
		The cuvette was washed once with water and dried by shaking off most of the water		
	12:32	From here on AE-441-2 and work under air. The flask was kept closed, when not in use	-2 before takeing sample.jpg	no change
		0.5 mL of AE-441-2 were transferred into the cuvette using a 1 mL Eppendorf syringe		
		To the cuvette water (1.5 mL) was added using a 1 mL Eppendorf syringe		
	12:35	An UV/Vis was recorded AE-441-2-1	-2-1.jpg	Looks ok, bit to concentrated
		The sample was pured in a 10 mL snap on cap vial		
		The cuvette was washed once with water and dried by shaking off most of the water		
		0.3 mL of AE-441-2 were transferred into the cuvette using a 1 mL Eppendorf syringe		
		To the cuvette water (1.8 mL) was added using a 1 mL Eppendorf syringe		
	12:38	An UV/Vis was recorded AE-441-2-2	-2-2.jpg	Good intensity

Analysis

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Processed Data	Comparative Data	Used for evalutaion	Interpretation
05.02	12:21	AE-441-1-1	UV/Vis	Equipment - Avantes UV/vis spectrophotometer	H2O	AE-441-1-1.ABS8	AE-441-1-1_7420287SP.TXT AE-441-1-1.png	KRA-092: Radiation of [Ru(CO)(OH)2(PNPtBu)], 365 nm and white	UV_Vis_plotting_a.py	A[365 nm]: 2.645 Stong absorbtion, no clear bands observed
05.02	12:25	AE-441-1-2	UV/Vis	Equipment - Avantes UV/vis spectrophotometer	H2O	Ae-441-1-2.ABS8	Ae-441-1-2_7420287SP.TXT AE-441-1-2.png	KRA-092: Radiation of [Ru(CO)(OH)2(PNPtBu)], 365 nm and white	UV_Vis_plotting_a.py	A[365 nm]: 0.741 Stong absorprtion, 2 weak bands observed. For detailed discussion see -2-2

05.02	12:28	AE-441-1-3	UV/Vis	Equipment - Avantes UV/vis spectrophotometer	H2O	AE-441-1-3.ABS8	AE-441-1-3_7420287SP.TXT AE-441-1-3.png	KRA-092: Radiation of [Ru(CO)(OH)2(PNPtBu)], 365 nm and white	UV_Vis_plotting_a.py	A[365 nm]: 0.261 1 peak and 2 bands observed. For detailed discussion see -2-2
05.02	12:35	AE-441-2-1	UV/Vis	Equipment - Avantes UV/vis spectrophotometer	H2O	AE-441-2-1.ABS8	AE-441-2-1_7420287SP.TXT AE-441-2-1.png	KRA-092: Radiation of [Ru(CO)(OH)2(PNPtBu)], 365 nm and white	UV_Vis_plotting_a.py	A[365 nm]: 0.504 1 peak and 2 bands observed. For detailed discussion see -2-2
05.02	12:38	AE-441-2-2	UV/Vis	Equipment - Avantes UV/vis spectrophotometer	H2O	AE-441-2-2.ABS8	AE-441-2-2_7420287SP.TXT AE-441-2-2.png	KRA-092: Radiation of [Ru(CO)(OH)2(PNPtBu)], 365 nm and white	UV_Vis_plotting_a.py	A[365 nm]: 0.315 1 peak at 229 nm (A = 2.269). 1 band at 275 nm (A = 0.916) 1band at approx. 365 nm (A = 0.315)

Results

Result	File
Absorption coefficient: 2110 +- 100 M*cm ⁻¹	AE-441_abostpiron_coefficientsaturated_conc.xlsx
Saturation concentration: 1.44 +- 0.04 mM	AE-441_abostpiron_coefficientsaturated_conc.xlsx

Linked experiments

- [KRA-092: Radiation of \[Ru\(CO\)\(OH\)2\(PNPtBu\)\], 365 nm and white](#)

PBDF - [NB-JSC-143: Testing of BFD2 \(MZ-124-6\) DMSO and DMSO-H2O solutions](#)

Prep work - [AE-442: Degassing of Milli-Q water](#)

Two-photon - [AE-382: Synthesis of \[Ru\(CO\)\(OH\)2\(PNPtBu\)\] using Ag2O, 5 h](#)

Linked resource

Equipment - [Avantes UV/vis spectrophotometer](#)


Attached files

[AE-441_abostpiron_coefficientsaturated_conc.xlsx](#)

sha256: a4a2d1d14b89ff52091eeca41266ddbb0b4988738b3f1d6e5e13f0a922e93fff


[AE-441-1-2.png](#)

sha256: e34c3781aedcd8b8642733af28e3739024de5f9b1cf65707831494b3d2b17a0b




AE-441-1-3.png

sha256: 93200439bdf6d8d397647325c0b7ac46d6512e63496930694db037b972a3b5d9




AE-441-2-1.png

sha256: 2a97443b4a8db2420a6e6f4f37a6c7296cd7f4064a731031b7a9a05ecfd13a03




AE-441-1-1.png

sha256: b935883cc24f2bd085cec2b1b310bea05e558cd50b25fd72d2e2f4c79c65f720



AE-441-2-2.png

sha256: 42980c91c3a8efae4f7303258cb6ce04412577bf9a623f7f4f46523761ac5a32



UV_Vis_plotting_a.py

sha256: 5025ed73591afe9e04388bb8e167d4e883e9642940d6519dc2305502cd77da14


water.jpg

sha256: 00493c1bd6c7ab9ff60f311cf518f0f86ed97e4dddc4736f6989799a91c7b540



1-2.jpg

sha256: 4bfe67e704273e5b6a2122aaabc74e3b098ee1e76a8ef01112990c4b9b7c8858



left-2-right-1.jpg

sha256: 351245fc435ab62889a7233b8479fc3b9e8cd2a5bf1aaf577949758084702d17



1-1.jpg

sha256: 98a79070c923a2142aaad8687bcabd1f4eb6dad0e7261bb23bc2d3149311b3a2



2-before-taking-sample.jpg

sha256: a8f8ffc159d805756861899cd4dad8419dd39f2c5c5583cc7cc09a2a333df7f7



2-1.jpg

sha256: 6eb2a4aff2f280e546594c94cd8136eef8926fff5e2f73a577d17916e6949a87



2-2.jpg

sha256: 3474d2ab782ff32a72a5f5ea281555801ab725321761e9a6aed2db73b9541459



AE-441-2-2_7420287SP.TXT

sha256: 2bdd74108f8278149530ffaa3ec08e509866b2060304b51d947ab341003e587e

AE-441-2-2.ABS8

sha256: ca7e4fe3c36e912c9e09345e01e3a5ded46686fd43bd0a7f6a70d24c7ada9626

AE-441-2-1_7420287SP.TXT

sha256: 944badf024f8d4ff2e70303df99a2bf8113900ff1a9fe01b45a9072467c11cc6

AE-441-2-1.ABS8

sha256: 9c575079a63d5c795c090ce0029f4205c4084ee6db3701b575faae9189af5c5b

AE-441-1-1.ABS8

sha256: ea7b7b5fbb7404e4e6bf6ee7c06327fe915264e5a9c93a0efcd459cdc387bebb

AE-441-1-3_7420287SP.TXT

sha256: 0366a50d03b8876bb9beae27a0804e320f3e78d553d1177e7e84eb0292a24b14

AE-441-1-3.ABS8

sha256: 32c25277ea29e137c1ab942398cfdba6e310e54699cf8c1eba61519f1681c3ff

Ae-441-1-2_7420287SP.TXT

sha256: 37e60a17b57f15c4fe6263ccc63ad712985d3c3641997456daa5d34970fcc2be

Ae-441-1-2.ABS8

sha256: 609e30387ff9c4a3e358a74e4a005327007877acef5dd7e2ce4358c0702a2520

AE-441-1-1_7420287SP.TXT

sha256: f2a8fdea17dec0558e9a57b52a5b644b9c713123198c08ff9b10dedfeee854fe



Unique eLabID: 20250205-abfeea3185aff5d7554811f3afa6144d09c46eb9

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1707>