


AE-450: Synthesis of Ag₂(17)O from AgNO₃ and H₂(17)O

Date: 2025-02-18
Tags: Ag₂O AE reference procedure AGNO₃ 17O
Category: Two-photon
Status: Done
Created by: Alexander Eith

Reaction scheme/sample structure

Chemical Formula: HNaO
Molecular Weight: 39,99677



Chemical Formula: AgNO₃
Molecular Weight: 169,87220

Literature/reference experiments

Literature	Preperative_Inorganic_Chemistry_Volume_2_Brauer, Link to book, p.1037 and 10.15227/orgsyn.036.0046
Reproduction	/
Similar experiments	Two-photon - AE-443: Synthesis of Ag₂O from AgNO₃

Reagents

[Ag217O-syn.xlsx](#) for calcualtion

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [μl]
Silvernitrate	7761-88-8	/	2.26	1	384	384.15	169.872	/	/
Sodium hydroxide	1310-73-2	/	2.26	1	90.04	90.07	39.9971	/	/
H ₂ (17)O 10 % in H ₂ (16)O	/	/	/	/	/	/	18.20	0.998	155 + 95 + 5

Work-up and Analytical Reagents

Name	CAS Number / Experiment Number	Inventory number	Mass _{exp} [g]	Volume [ml]	Concentration [M]
Milli-Q water	/	/	/	approx. 70	/

Procedure/observations

Date	Time	Step	Pictures	Observations
18.02	13:35	In a GC vials with screw caps AgNO ₃ (flask 1) was weighed		both white solids
		In a GC vials with 300 µL micro inserts with screw caps NaOH (flask 2) was weighed	left NaOH right AgNO ₃ .jpg	
	13:50	To the flask 1 H ₂ (17)O 10 % in H ₂ (16)O (155 µL) and to the flask 2 H ₂ (17)O 10 % in H ₂ (16)O (95 µL) was added using a 100 µL Eppendorf pipette		
	13:54	flask 2 was vortexed (Equipment - VWR® VV3, Vortex Mixer , stage 5/6, than 6/6) for approx. 3 min and than approx. 4 min		some solids remain at the bottom
	13:55	flask 1 was vortexed (Equipment - VWR® VV3, Vortex Mixer , stage 5/6) for approx. 1 min	left NaOH right AgNO ₃ after addition of water.jpg	clear solution
	14:05	to flask 2 H ₂ (17)O 10 % in H ₂ (16)O (5 µL) was added using a 10 µL Eppendorf pipette		
	14:04	flask 2 was vortexed (Equipment - VWR® VV3, Vortex Mixer , than 6/6) for approx. 1 min	after vortexing left NaOH right AgNO ₃ .jpg	some solids remain at the bottom
	14:06 - 14:09	The mixture in flask 2 was taken up in a 100 µL Eppendorf pipette for approx. 30 times to yield dissolution of remaining solids	NaOH after Eppendorf.jpg	slightly cloudy solid
	14:09	The solution of flask 2 was added to flask 1 using a 100 µL eppendorf pipette	After addition of NaOH.jpg	after addition brown solid forms at layer
	14:12	The obtained mixture was vortexed (Equipment - VWR® VV3, Vortex Mixer , stage 2/6) for approx. 20 and than at stage 4/6 for approx. 2 min. In between the sample was shaken (approx. 3 times for 5 s)	after vortexing.jpg	after vortexing more brown solid forms
		The vial was stored at rt in the dark	after reaction.jpg	bit darker, more blackish, mixture solidified
	15:35	To the mixture Milli-Q water (approx. 1 mL) was added using a squirt bottle		
		The mixture was vortexed (Equipment - VWR® VV3, Vortex Mixer , stage 4/6) for 10 s		dark suspension

		The suspension was transferred into a 14 mL centrifuge tube using a glas pipette		
		This was repeated with the 1 mL and 1.5 mL of water		
	- 15:40	The centrifuge tube was filled to 10 mL with Milli-Q water		
	15:45 -	The mixture was centrifuged (Equipment - Centrifuge Hettich EBA 20 , 2 min, 6000 rpm) and the supernatent liquid was taken off using a glas pipette	after first centrifugation.jpg	clear liquid, dark solid
	- 16:15	The obtained solid was washed with Milli-Q water (6 * 10 mL) by centrifuged (Equipment - Centrifuge Hettich EBA 20 , 2 min, 6000 rpm) and taken off the supernatent liquid	combined liquid phases.jpg	no changes in appreance
	16:20	After the last centifugation water (approx. 1 mL) was added in the tube and the suspension was transffered into a 5 mL snap on cap vial using a glas pipette. To get a good suspension the liquid was taken up into the pipette several times		
		To the tube water (0.3 and 0.1 mL) was added and the obtianed suspension was transffered into the vial		
		The supernatant soultion in the vial was taken off into the tube and shortly centrifuged.		
		Ther supernatant soultion in the tube was decanted and the remaining suspension was transffered into the vial	solids.jpg	balck slurry
	16:35	The top of the vial was covered with puntucered aluminium foil and placed in an drying oven over night at 95 °C		
19.02	07:50	The vial was removed from the oven capped and stored at room temerpature in the dark.		
20.02	09:00 - 12:40	The solid was dried under reduced pressure (Schlenk line) at rt		
		The solid was stored in the dark at rt		

Analysis

Product characterization

Sample	Mass [mg]	Purity	Mass _{pure} [mg]	Amount [mmol]	Yield [%]	Description	Image	Storage location
AE-450	225.02	/	/	0.97	86	dark dray, black solid	product.jpg	E004 SSC

Results

Product obtained in good yield

Linked experiment

Two-photon - [AE-443: Synthesis of Ag₂O from AgNO₃](#)

Linked resources

Equipment - [VWR® VV3, Vortex Mixer](#)

Equipment - [Centrifuge Hettich EBA 20](#)

Attached files

AE-443.cdxml

sha256: 9fd79b07a2c2244119f5827b96960f021b248a9d4b21eaea37527e0d6b109da6

Ag2170-syn.xlsx

sha256: 3c7ec664399b84e1dc5b5efc5507483fa3eeb8d60b01513db206cc764f0ea286

after-reaction.jpg

sha256: df77c0f0ac46fa7a634326b79d5115e39a9bd70914232bf5ace4b6ec69765617



after-first-centrifugation.jpg

sha256: 193579c75eb3643e43599949583bb8d264c20027bcb83ade0c243a5892c8466d



solids.jpg

sha256: d748655b50a43e2d15dc8e64d480f4299c1f1038c7ebdc53b0413f104273300c



product.jpg

sha256: 9dfcc7ab7b4658499edc041769efaa8d729e042790146f31e709c16472fa1d6c



combined-liquid-phases.jpg

sha256: bf2d079873af466715282c6d50200dc7f0a44dfd6e8d853db59eb54a8f7b0566



after-vortexing.jpg

sha256: 87435274efe2c77861a69ad47e7bb1f4cbc53a55a47ab7c94e0074e5ce5e1a7d



left-NaOH-right-AgNO3.jpg

sha256: fa6f073ce84fa36bfae1afbc96056d3c5e6cfe99ceb021535943aa238387e801



left-NaOH-right-AgNO3-after-addition-of-water.jpg

sha256: 9f8b371550bb559fedd2d6c49124250daf3c0d0ff49982e3e0b375627a3b6f7e



NaOH-after-Eppendorf.jpg

sha256: e8b56a30b868f13d80a9b634079f2b33e5a0490d366d8fe68783d2518fb0fbe6




after-vortexing-left-NaOH-right-AgNO3.jpg

sha256: d99054b5896cffe45577c95954ab0b836adbf7e256b488b7463ffb4467a65612



After-addition-of-NaOH.jpg

sha256: e6538de76f568384008b5fd8b9794e110b1f0a3fb942dc02d4bb79ecc24aa4f9



Unique eLabID: 20250218-1466710a3583ee1cc1a11ee0b65f2504a4dfd2d8
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1771>