

AE-592: Dual irradiation of [Ru(CO)(OH)₂(PNP)] (AE-570)

Date: 2025-10-14

Tags: Dihydroxo PNP P(tBu)N(py)P(tBu)
[Ru(CO)(OH)₂(PNP)] NMR AE Future ¹H
³¹P Dual irradiation two-photon
photocatalysis

Category: Two-photon

Status: Done

Created by: Alexander Eith


Objective

Testing dual irradiation effect of [Ru(CO)(OH)₂(PNP)] at higher irradiance, to get better conversion after 16 h

Reproduction of AE-590

Reaction scheme/sample structure

Radiation of Dihydroxo



ChemDraw File (linked): [AE-590.cdxml](#)

Literature/reference experiments

Literature	Procedure after https://doi.org/10.1039/D1EE01053K
Reproduction	Two-photon - AE-590: Dual irradiation of [Ru(CO)(OH) ₂ (PNP)] (AE-570)
Similar experiments	Two-photon - AE-519: Dual Radiation of [Ru(CO)(OH) ₂ (PNP{tBu})] (AE-456-1) Two-photon - AE-511: Dual Radiation of [Ru(CO)(OH) ₂ (PNP{tBu})] (AE-456-1) - reproduction II Two-photon - AE-502: Dual Radiation of [Ru(CO)(OH) ₂ (PNP{tBu})] (AE-456-1) - reproduction Two-photon - AE-500: Dual Radiation of [Ru(CO)(OH) ₂ (PNP{tBu})] (AE-456-1)

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]
------	--------------------------------	------------------	---------------	-------------	---------------------------	--------------------------	--------------------	----------------	-------------	--------------------

[Ru(CO)(OH)2(PNP{tBu})]	Two-photon - AE-591: Stocksolution of [Ru(CO)(OH)2(PNP{tBu})] in water	/	/	/	/	/	558.65	/	3 * 0.6	concentrated (approx. 1.44 mM)
-------------------------	---	---	---	---	---	---	--------	---	---------	--------------------------------------

Irradiation Parameters

Power measurement was performed in experiment [Two-photon - AE-590: Dual irradiation of \[Ru\(CO\)\(OH\)2\(PNP\)\] \(AE-570\)](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation set-up number	Equipment - Irradiation setup 2 (CEEC II, E004)

	Light Source Name	Used power source	Wavelength [nm]	Power Setting [mW]	Analogue Setting [-]
First light source	Light Source - UHP LED 365 nm-2	Power Sources - BLS-18000-1 1	365	500	9.68
Second light source	Light Source - UHP LED 470 nm	Power Sources - BLS-13000-1E	470	600	8.94

AE-592-1

Used beam combiner [Name or None]	Beam Combiner - Beam Combiner LCS-BC25-0409 -1
Used light source	365 + 470 nm
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	25
Stirring speed [rpm]	/
Irradiation duration [h:min]	17:32

AE-592-2

Used beam combiner [Name or None]	Beam Combiner - Beam Combiner LCS-BC25-0409 -1
Used light source	365 nm
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	25
Stirring speed [rpm]	/
Irradiation duration [h:min]	17:30

AE-592-3

Used beam combiner [Name or None]	Beam Combiner - Beam Combiner LCS-BC25-0409 -1
Used light source	470 nm
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	25
Stirring speed [rpm]	/
Irradiation duration [h:min]	17:28

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard [Protocol - Schlenk Technique](#)

The synthesis/reactions were performed in the dark by either wrapping the reaction flasks in aluminum foil or turning the lab light and light in the fume hood off, when the flasks were uncovered.

Date	Time	Step	Observations	Pictures
13.10	16:05	Two Young type NMR tube was prepared according to Protocol - Preparation of NMR Sample (from step 3) using 0.65 mL of AE-591	AE-592-1 and AE-592-2	/

		From AE-592-1 an NMR was measured (AE-592-1-1)	AE-592-1-1 , slight yellow soultion, maybe very small amount of solid deposited on the wall	AE-592-1-1.jpg
		From AE-592-2 an NMR was measured (AE-592-2-1)	AE-592-2-1 , slight yellow soultion	AE-592-2-1.jpg
14.10	17:14	The power was adjusted using Power Meter - 843-R-USB + 919P-020-12 in Equipment - Advanced power measurment setup V1.0 (365 nm: 500 mW, 470 nm: 500 mW)	/	/
	17:19	AE-590-1 was placed in the setup using the Photoreactor - AE-558: Design and print of holder and holder disk for 2 NMR tubes Holder_2NMR_tubes	no change	back before irrad 1.jpg
	17:20	The irradiaiton of AE-592-1 was started	/	back during irrad 1.jpg
	10:52	The irradiaiton of AE-592-1 was stopped	/	back after irrad 1.jpg
	12:42	From AE-592-1 an NMR was measured (AE-592-1-2)	AE-592-1-2 , darker orange, some solids formed on the walls	AE-592-1-2.jpg
	14:20	One Young type NMR tube was prepared according to Protocol - Preparation of NMR Sample (from step 3) using 0.65 mL of AE-591 using a Prep work - AE-489: Capillary production	AE-592-3	/
	14:50	From AE-592-3 an NMR was measured (AE-592-3-1)	AE-592-3-1 , slight yellow soultion	AE-592-3-1.jpg
	16:41	AE-592-2 was placed in the setup using the Photoreactor - AE-558: Design and print of holder and holder disk for 2 NMR tubes Holder_2NMR_tubes	no change	back before irrad 2.jpg
	16:42	The irradiaiton of AE-592-2 was started	/	back during irrad 2.jpg
	10:12	The irradiaiton of AE-592-2 was stopped	no significant change	back after irrad 2.jpg
	10:43	From AE-592-2 an NMR was measured (AE-592-2-2)	AE-592-2-2 , no change	AE-592-2-2.jpg
	15:05	AE-592-3 was placed in the setup using the Photoreactor - AE-558: Design and print of holder and holder disk for 2 NMR tubes Holder_2NMR_tubes	no change	back before irrad 3.jpg

	15:10	The irradiation of AE-592-3 was started	/	back during irradiation 3.jpg
	08:38	The irradiation of AE-592-3 was stopped	no significant change	back after irradiation 3.jpg
	09:53	From AE-592-3 an NMR was measured (AE-592-3-2)	AE-592-3-2 , no change	AE-592-3-2.jpg

Analysis

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
14.10	11:22	AE-592-1-1	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H ₂ O	AE-592-1-1.zip	/	AE-592-1-1_10.nmrium	Two-photon - AE-590: Dual irradiation of [Ru(CO)(OH) ₂ (PNP)] (AE-570)	As expected, before irradiation, reactant, and minor impurity
15.10	16:27	AE-592-1-2	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H ₂ O	AE-592-1-2.zip	/	AE-592-1-1_10.nmrium	Two-photon - AE-590: Dual irradiation of [Ru(CO)(OH) ₂ (PNP)] (AE-570)	Mainly as expected: formation of hydride species and side products, but more side products as expected
14.10	15:47	AE-592-2-1	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H ₂ O	AE-592-2-1.zip	/	AE-592-1-1_10.nmrium	/	As expected, before irradiation, reactant, and minor impurity
16.10	14:17	AE-592-2-2	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H ₂ O	AE-592-2-2.zip	/	AE-592-1-1_10.nmrium	/	No conversion observed
16.10	02:02	AE-592-3-1	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H ₂ O	AE-592-3-1.zip	/	AE-592-1-1_10.nmrium	Two-photon - AE-500: Dual Radiation of [Ru(CO)(OH) ₂ (PNP{tBu})] (AE-456-1)	As expected, before irradiation, reactant, and minor impurity
17.10	14:32	AE-592-3-2	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H ₂ O	AE-592-3-2.zip	/	AE-592-1-1_10.nmrium	Two-photon - AE-500: Dual Radiation of [Ru(CO)(OH) ₂ (PNP{tBu})] (AE-456-1)	As expected, no conversion observed

Results

Sample	Yield (hydride species, cis) [%]	Yield (hydride species, trans) [%]	Yield (dihydroxo) [%]	Yield (main side product, 57 ppm) [%]	Yield (soluble side products) [%]	Yield (other side products) [%]	Data
AE-592-1	17	0	7	41	20	15	AE-592.xlsx
AE-592-3	0	0	88	0	3	9	AE-592.xlsx

AE-592-2: no evaluated, since LED was connected to wrong port of the beam combiner

Worked in general, good agreement of AE-592-1 with AE-590-1 for hydride yield, but more sideproducts formed

AE-592-3: as expected no hydride formation

Future recommendations

Old procedure	Problem	Suggested new procedure
LEDs connected to the wrong ports of the beam combiner for irradiation with 365 nm	way less power output than desired	do connect them correctly

Linked experiments

Two-photon - [AE-500: Dual Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\]](#) (AE-456-1)

Two-photon - [AE-502: Dual Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\]](#) (AE-456-1) - reproduction

Two-photon - [AE-511: Dual Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\]](#) (AE-456-1) - reproduction II

Two-photon - [AE-519: Dual Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\]](#) (AE-456-1)

Two-photon - [AE-581: Stock solution of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\]](#) in water

Two-photon - [AE-590: Dual irradiation of \[Ru\(CO\)\(OH\)2\(PNP\)\]](#) (AE-570)

Two-photon - [AE-591: Stock solution of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\]](#) in water

Linked resources

Beam Combiner - [Beam Combiner LCS-BC25-0409 -1](#)

Equipment - [Advanced irradiation setup V1.0 I](#)

Equipment - [Advanced power measurment setup V1.0 I](#)

Equipment - [Manual irradiation setup](#)

Equipment - [Irradiation setup 2 \(CEEC II, E004\)](#)

Light Source - [UHP LED 470 nm-1](#)

Light Source - [UHP LED 365 nm-2](#)

Power Meter - [843-R-USB + 919P-020-12](#)

Power Sources - [BLS-13000-1E](#)

Power Sources - [BLS-18000-1 1](#)

Attached files

AE-592.xlsx

sha256: 7a75d51fa97d08e3c0e9bc4e82021cdd7e93699eb5b9439527dd6cdabab72f37

AE-592-1-1_10.nmrium

sha256: 4801770f68e900bcfa702276485d50cefa7c1579455f2e461ba9a03d01dbb799

AE-592-3-2.zip

sha256: 992a366b7ecd5262b1ee8acc0517c0902597b5d26d24b7ff2022feb1e948b

AE-592-3-1.zip

sha256: dcfffd7054d3fe265aba04af94d1fc9026133bc8d34b2b9426e17bf8fa02dd3b

AE-592-2-2.zip

sha256: dc92855820457aff12b07656c31f4db1e50966b4a53dd8e014da30c70e6ba79f

AE-592-2-1.zip

sha256: a96f27d5e644ea07a29a17a1a5bdea0ddcfba8e65b4b0f6f5c225e6d1aff302c

AE-592-1-2.zip


sha256: 137154fab1e4f513f7f94cd5be0540d35733e348821dc1d72d8aa8c5654dfb82

AE-592-1-1.zip

sha256: fa43f5c3df31dbc1451208857f27ab21f2c10ecf31a9bc69d9c3a01e4f06a337

AE-590.png

sha256: 4acf7b189a4440e933f6106b5439aec65f140eb6486dd99d23bb59f337309365



AE-590.cdxml

sha256: 2bb17a23d6ac9133244c0395480ff7b2d05f24ff8c1a41893e7170ba2481ccdc

AE-592-3-2.jpg

sha256: c1a54b03a9a33eeeedec63141ce821d90554c3eca56cdd3ca1c4b8da62dc8afb



back after irrad 3.jpg

sha256: eab1cebbde97421b8e3b51f76b3d86de0cc7f208c2bf22f7a3f546abb4f35fea



back during irrad 3.jpg

sha256: 49ef822f4fd18fc7ef096086b06d427d525b63e05f3b13e15421f7adf69b674a



back before irrad 3.jpg

sha256: 499435fd030128725c45016814d781a837db09ffae93190ac46d3a122519537



AE-592-2-2.jpg

sha256: b4796607435e9170abe06a99b3ff8c3dfacbd8962ae7fcf038e9a96f7fa8a41c



back after irrad 2.jpg

sha256: 4c2ba3cdad1a05cbd64714c0a9fc238ddb5a7e5ccd7ea2d8285ae5bea73e90f3



back during irrad 2.jpg

sha256: 57609f6d810cc5ce15faa0bac092ca5d2db698df891724ec59b1bb6240970201




back before irrad 2.jpg

sha256: 1c3d874d3273610fa0b744ff1641a51dfb48e28fef3e3623af22fb22f5e4aa46




AE-592-3-1.jpg

sha256: d3fc3320fe6718b6bf88782cfd0461f8b3ac5532407709ca586964cf53a27418



AE-592-1-2.jpg

sha256: da63bd6b5ed8069bf7e468e0db870e91ba2a76304d1e34d30223226754295ca6



back after irrad 1.jpg

sha256: 3343a14ae5c87b9b9f8b7abbbff8a5fd965343e3c607ff7011726e654eb7fdaa2



back during irrad 1.jpg

sha256: e390c2d34cda268f485f1fb15c1ab7e3da539c227b1dc7c42fdcf1d1e4aad47



back before irrad 1.jpg

sha256: c9f72ac851e0db3c78df01001b031de6aadd9bf4dde531134f820ac78d62ebd7




AE-592-2-1.jpg

sha256: d08794e6a0f0383becdb735f5f23e6846d4d2b77d0c57b86bdf75d73bf31207b



AE-592-1-1.jpg

sha256: a0a4def417c5eae492833542c810163790f5acd335a07d36616202cf7d717203



AE-590.cdxml

sha256: 7110161878b2b395970d517d41146467e0d34d6aba8df3850834b3a575167a09



Unique eLabID: 20251014-03590b6be1b3c11751aef73b6c0273b4380e8c67

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3184>