

Irradiations using Young type NMR tube

Date: 2024-09-13
Tags: NMR AE Irradiation
Category: Protocol
Created by: Alexander Eith

Goal

Irradiation of Young type NMR (typically for catalytic studies with NMR)

Prerequisites and preparation

- Young type NMR tube

Equipment - Irradiation Set-Up

Note: before doing the experiment thing about apropiate light sources (best measure UV/Vis beforehand of compound), irradiation intensity (sun has approx. 100 mW/cm², sensor has area of 1.13 cm²) and duration

Steps

Step number	Step description
1.	Prepare NMR sample (use Protocol - Preparation of capillary for H2O-NMR as external standart if possible, for exact quantification an external standart is needed)
2.	Measure NMR before irradiation
3.	Build irradiation setup according to Protocol - Irradiation setup
4.	Place NMR tube in setup and start irradiation
5.	After irradiation measure an NMR

Linked resources

Archived - [Irradiation Set-Up](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 1\)](#)

Protocol - [Preparation of capillary for H2O-NMR](#)

Protocol - [Irradiation setup](#)

Protocol - [O2 Measurement procedure](#)



Unique eLabID: 20240913-26f186bab040c3aa32e5c0820b15d06d2b87271d

Link: <https://elab.water-splitting.org/database.php?mode=view&id=140>