

Goal

This protocol should highlight the important tricks, when using J. Young solvent flask.

There are two different sizes: 500 ml (big ones) are used for dry solvents and the 25 ml (small ones) are used for dry NMR solvents.

Filling the J. Young flask


If the solvent is taken out of the SPS: [Protocol - Collecting of SPS solvents](#)

Degassing, when using this type of flasks: [Protocol - Degassing of organic solvents in Young-type flask](#)

Setup

First steps to take, when using J. Young type flasks.

Step number	Step description
1	Place a septum on the opening
2	Place tissues under the flask
3	Clamp the flask under the "Y-thing"
4	Connect the flask to the line using GL14



Taking solvent out of the J. Young flask using a syringe

Here, the way without touching/damaging the PTFE valve is described. Important: the tip needle should not touch the PTFE valve!

There is a video of the whole process: [whole_process_01-05.mov](#)

And two videos of the steps 06-07, where the movement of how NOT to touch the PTFE valve is shown precisely: [steps_04_05.mov](#), [steps_04_05_other_angle.mov](#)

Step number	Step description
Preparation	See "set up" above setup.jpeg and picture preparation.jpeg
01	Preparation of the needle, the needle is bend at the top part (ca. 1 cm from the end) in ca. 35 °. For this always use a tissues, since you don't know what is on your gloves. 01_needle.jpeg , 01_needle_closer.jpeg

02	After three cycles of changing the atmosphere in the needle, place the needle through the septum. 02_going_in.jpeg Go with the needle up to the PTFE valve, but still keep a distance. <i>Important: the tip needle should not touch the PTFE valve!</i> 02_1_first_stop.jpeg
03	Open the PTFE valve up to the point where it is fully open, but the ring (bottom part) stays on. 03_opening.jpeg
04	Go through the opening next the the PTFE valve in a way that the tip of the needle is sliding above the glas and pointed downwards. <i>(Better seen in the video)</i> 04_going_through.jpeg Video: steps_04_05.mov, steps_04_05_other_angle.mov
05	Now just push the needle as much as you can. 05_you_are_in.jpeg Video: steps_04_05.mov, steps_04_05_other_angle.mov
06	Turn the needle by 180 ° and you should be able to take of the solvent. For this the flask needs to be clamped out and tilted as much as you need to. 06_taking_the_solvent.jpeg <i>Since these are solvents, you should omit putting the solution back to the solvent flask and just take everything out what was once in your needle.</i>
07	Take the needle above the "Y"-part and bend the needle using your fingers. <i>Never bend needles using the septum.</i>
08	Close the PTFE valve.
09	Take an argon puffer on top of your solution in the needle and take it out.
10	Congratulation, you are done :)

Linked resources

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [Collecting of SPS solvents](#)

Attached files

setup.jpeg

sha256: ed9a79e933167845fdccb1a98aa7c124c6c00d3ebc6e6f971cabeb4a89e9f976



whole_process_01-05.mov

steps_04_05.mov

sha256: 1397de03d0249a013f9a0730e1f479650c2d667f36e82df5d3b03bc910bd7e2a

steps_04_05_other_angle.mov

sha256: fc9774783cd12a7cb06c1865e78bde3519986384c3dbc3a51244b472565a3225

04_going_through.jpeg

sha256: 2d3741ac092dd12293df557fd2867c80619fef6e08c69b5cad6d1e2daa2d4810



05_you_are_in.jpeg

sha256: bdacc327b52f80967bdb817278251c49f552b826126efcc867f5d0fb9c6784f9



06_taking_the_solvent.jpeg

sha256: 3d490a132fa2e1b31b290e8ec3dbaca4199cebcc9d0a2929e1d3255e8d06eb6d



01_needle_closer.jpeg

sha256: 780e8691254d7f8c7419724ef8eb7b9768f30910b7e919a2dd4dc2bcd7286cf4



02_going_in.jpeg

sha256: f92c3e8e99bb92728baf68aadbab95d1fa16b329f4840a5dca611a29609a2218



02_1_first_stop.jpeg

sha256: 6e1c45031df251635070f538e3a1e89d82dca1769f57c10a0be47d3ad9d6505f



03_opening.jpeg

sha256: 835b3aa0c07d9d19bc97b69c4d379399a835affdf21b11aac38ae0a59266df3c



preparation.jpeg

sha256: 67c9c79f7fd81e3fa521bf2669918038f0c480b95bbc08b64bd560403127fbad



01_needle.jpeg

sha256: eccef27eb96cbeeedecbaf732ef836f6aaa078f909c592af13e9a1e3b7a08a32



Unique eLabID: 20240909-6901bb124449cf2ba233bfa23d7508c34d40cd84
Link: <https://elab.water-splitting.org/database.php?mode=view&id=132>