

AE-120: Degassing of n-Heptane

Date: 2023-08-02

Tags: Degassed AE

Status: Done

Created by: Alexander Eith

Reaction scheme/sample structure

Reagents

Name	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
n-Heptane (Reactpur, VWR, 99 %)	/	/	/	/	/	approx. 450 mL

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard Schlenk technique or in an argon filled glovebox (MBRAUN).

Date	Time	Step	Observations
02.08	10:50	A 500 mL-solvent flask was prepared.	
		Molar sieve 3 A (approx. 100 mL) was filled in the flask and the atmosphere was subsequently changed to argon.	
24.08	8:35 - 8:50	n-Heptane was transferred into the flask according to [Protocol] Cannula Filtration/Transfer	
	9:05 - 9:45	Argon was sparged through the solution analog to [Protocol] Degassing of fluids for 40 min	

Analysis

Date	Time	Sample name	Analysis method	Solvent	Raw Data	Processed Data	Interpretation

Product characterization

Linked resources

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [KFT \(Karl-Fischer-Titration\)](#)

Protocol - [Collecting of SPS solvents](#)



Unique eLabID: 20230731-3c64d105288bb4e528901d60723ed1a068f99fd1
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=567>