

# AE-382: Synthesis of [Ru(CO)(OH)<sub>2</sub>(PNP*t*Bu)] using Ag<sub>2</sub>O, 5 h

Date: 2024-11-13


Tags: Dihydroxo Ag<sub>2</sub>O PNP  
P(*t*Bu)N(*py*)P(*t*Bu) [Ru(CO)(OH)<sub>2</sub>(PNP)]  
NMR ESI-MS AE MS reference analytics  
1H 31P Elemental Analysis

Category: Two-photon

Status: Done

Created by: Alexander Eith

## Reaction scheme/sample structure



## Literature/reference experiments

Literature	<a href="https://doi.org/10.1039/D1EE01053K">https://doi.org/10.1039/D1EE01053K</a>
Reproduction	/
Related experiments	Experiment - AE-KRA-204: Synthesis of [Ru(CO)(OH) <sub>2</sub> (PNP <i>t</i> Bu)] using Ag <sub>2</sub> O, 5 h Experiment - AE-238: Synthesis of [Ru(CO)(OH) <sub>2</sub> (PNP <i>t</i> Bu)] using Ag <sub>2</sub> O, 5 h Experiment - AE-252: Synthesis of [Ru(CO)(OH) <sub>2</sub> (PNP <i>t</i> Bu)] using Ag <sub>2</sub> O, 5 h AE-310: Synthesis of [Ru(CO)(OH) <sub>2</sub> (PNP <i>t</i> Bu)] using Ag <sub>2</sub> O, 5 h

## Reagents

Name	CAS-Number/Expierment-Number	Amount [mmol]	Equivalen ts	Mass <sub>theo</sub> [mg]	Mass <sub>exp</sub> [mg]	Molar mass [g/mol]	Volume <sub>theo</sub> [ml]	Volume <sub>exp</sub> [ml]	Density [g/mL]
[RuCl <sub>2</sub> (CO)(PNP <i>t</i> Bu)]	Experiment - AE-261: Synthesis of [RuCl <sub>2</sub> (CO)(PNP)] with 2,6-Bis(tert-butyl)phosphinomethylpyridine via [RuCl <sub>2</sub> (CO)(p-cymene)] (AE-261-1)	0.827	1.00	-	492.57	595.53	/	/	/
Ag <sub>2</sub> O	20667-12-3	5.078	6.14	1176.88	1176.70	231.74	/	/	/
H <sub>2</sub> O	Experiment - KRA-104: Degassing of Milli-Q water	193.54	234	3487	/	18.02	3.49	3.6	0.998
THF	Experiment - KRA-080: Degassing of THF	/	/	/	/	/	25	25	/
n-Heptane	Experiment - AE-121: Degassing of n-Heptane	/	/	/	/	/	/	20	/

## Work-up Reagents

Name	CAS Number /	Inventory number	Mass <sub>exp</sub> [g]	Volume [ml]	Concentration [M]
THF	Experiment - KRA-080: Degassing of THF	/	/	25 +20	/
n-heptane	Experiment - AE-121: Degassing of n-Heptane	/	/	59 + 55	/
MeCN	AE-152: Degassing of MeCN	/	/	6	/

## Procedure/observations

All steps, unless mentioned otherwise, were carried out under argon using standard Schlenk technique or in an argon filled glovebox. (see Protocol [Protocol - Schlenk Technique](#))

Due to the instability of Ag<sub>2</sub>O and Dihydroxo to the light, the procedure was done mostly in the dark, by either wrapping the flasks containing either or both in Aluminium foil or by tuning off the light in the lab and fumehood

Date	Time	Step	Observations
14.11	9:10	A 50mL- Schlenk flask was prepared, by heating under vacuum	
	9:30	[RuCl <sub>2</sub> (CO)(PNP)] and Ag <sub>2</sub> O was weighted in air into the reaction schlenk flask and the flask was set under argon.	
		The light in the fumehood was turned out	
	10:10	THF (25 mL) was added to the reaction schlenk flask.	
	10:15	Water (3.6 mL) was added to the reaction schlenk flask	
	10:15 - 15:15	The reaction mixture was stirred for 5 h at 700 rpm. During this the flask was wrapped with aluminium foil	<a href="#">Start of reaction.jpg</a> <a href="#">after reaction.jpg</a>
	15:15 - 15:20	The reaction mixture was allowed to settle	

	15:20 - 15:50	The mixture was filtered through a pore 4 frit according to <a href="#">Protocol - Filtration with frit technique</a> into a 100 mL-Schlenk	
	15:50 - 16:05	The reaction flask was washed with MeCN (2 * 3 mL)	1st washing: yellowish liquid phase 2nd washing: nearly clear liquid phase
	16:05 - 17:05	The solvents of the filtrate were removed under reduced pressure. • When approx. 5 mL were left the water bath was removed and the mixture was allowed to freeze (at approx. 16:20) • At 16:45 the water bath was added again to slowly melt the remaining solvent	This procedure avoids bubbling up of the mixture
		The solid was stored in the dark at - 20 °C	Dark brownish solid <a href="#">Solid after night.jpg</a>
15.11	10:20 - 10:30	To the solid THF (25 mL) was added and the obtained mixture was stirred and shaken strongly for approx. 7 min	<a href="#">after adding THF.jpg</a>
	10:30 - 10:37	n-heptane (50 mL) was added to the mixture	<a href="#">after adding n-hept.jpg</a>
	11:17 - 12:00	The mixture was filtered through a pore 3 frit according to <a href="#">Protocol - Filtration with frit technique</a> into a 100 mL-Schlenk	
	12:03 - 12:10	The solid was washed with n-heptane (3 * 3 mL)	
	12:13 - 13:17	The solid was dried under reduced pressure	
	13:30 - 13:40	The solid was transferred into a 2 mL vial in a 10 mL Schlenk flask	Some material was lost, since it stuck strongly to any surface
		During this an aliquote (approx. 5 mL) was transferred into another 10 mL Schlenk flask	
	13:45 - 14:35	Both flasks were dried under reduced pressure	
		Both flask were set under argon	
	14:50	From the aliquote an NMR sample was prepared in DMSO-d6 according to <a href="#">Protocol - Preparation of NMR Sample</a>	<b>AE-382-1</b> <a href="#">AE-382-1.jpg</a>
		The remaining solid was stored in the dark under argon at -20 °C	<b>AE-382-1</b> <a href="#">Product.jpg</a>

04.12		Further workup, since in <sup>1</sup> H NMR some more sideproducts are observed (with all remaining <b>AE-382-1</b> )	
		A 100 mL flask was prepared	
		The vial containing <b>AE-382-1</b> was added to the flask is Ar counterflow	
	16:05	To the 100 mL flask THF (20 mL) was added	<a href="#">in THF.jpg</a> dark suspension
		It was tried to remove th vial without sucess	
	16:07 - 16:12	The flask was shaken to get all solids out of the vial	
	16:12 - 16:23	The flask was stirred	brownish suspension
	16:23 - 16:30	To the flask n-heptane (52 mL) was added	<a href="#">after n-hept.jpg</a> orange solution, brownish.gray solid
	16:30 - 16:32	The suspension was stirred	
		The flask was stored in the dark in <a href="#">Equipment - Freezer Lab E004 CEEC II</a> under Ar	
05.12		The obtained suspension was warmed to rt	<a href="#">before filt.jpg</a> orange solution, brownish.gray solid
	10:40 - 14:30	The suspension was filtered through a pore 4 filter frit	Very slow filtration, filtration was stopped after approx. 40 mL of solvent
	14:30 - 14:50	The obtained filter residue was washed with n-Heptane (1 * 1.5 mL)	<a href="#">washing of filterate on 1st frit.jpg</a>
	14:50 - 14:58	The obtained solid ( <b>AE-382-2</b> ) was dried under reduced pressure	
	14:53	The filtration setup was changed to a pore 3 frit. Therefore the new frit was connected with a "Vakuumvorstoß", closed and set under inert conditons. Afterwards the cannula, still conected to the flask containing the suspesion, was attached in Ar counterflow. Then the reciving flask was added under Ar counterflow and the filtration was continued	
	14:56 - 15:13	The remaining suspension was filtered through a pore 3 filter frit	
	15:15 - 15:20	The obtained residue ( <b>AE-382-3</b> ) was washed with n-heptane (3 * 1 mL)	

	15:30 - 17:05	AE-382-2 and -3 were dried under reduced pressure	
		Both were stored on the frit in the dark under Ar in <a href="#">Equipment - Freezer Lab E004 CEEC II</a>	
09.12	10:40	<b>AE-382-3:</b> The solid was transferred into a 1.5 mL screw cap vial in a 10 mL Schlenk flask. during this approx. 5 mg were added to another 10 mL Schlenk flask to prepare an NMR sample.	
		Both flasks were set under argon and stored in <a href="#">Equipment - Freezer Lab E004 CEEC II</a> in the dark under Ar	
12.12	11:05	To the flask containing the NMR sample DMSO-d6 (0.75 mL) was added	<a href="#">-3 in DMSO-d6.jpg</a> brownish solution
	11:05 - 11:10	The obtained solution was stirred	
	11:10	<a href="#">An NMR sample was prepared according to Protocol - Preparation of NMR Sample</a>	<a href="#">AE-382-3.jpg</a> brownish solution
15.01	10:05	<b>AE-382-2:</b> The solid was transferred into a 1.5 mL screw cap vial in a 10 mL Schlenk flask. during this approx. 4 mg were added to another 10 mL Schlenk flask to prepare an NMR sample.	
		Both flasks were set under argon and <b>AE-382-2</b> was stored in <a href="#">Equipment - Freezer Lab E004 CEEC II</a> in the dark under Ar	
	10:25	To the flask containing the NMR sample MeCN-d3 (0.65 mL) was added	
		<a href="#">An NMR sample was prepared according to Protocol - Preparation of NMR Sample</a>	<a href="#">AE-382-2.jpg</a>
13.03		A 10 mL and a 5 mL Schlenk flask were prepared, each containing a 1.5 mL screw cap vial without lid	
	13:45	An sample for an elemental analysis was prepared by weighing in 3.75 mg of <b>AE-382-2</b> into the vial of the 10 mL Schlenk flask and setting it under argon afterwards.	<b>AE-382-EA</b> (AE382a on sheet, since max. 6 letters)
	13:50	An sample for MS was prepared. Therefore 1.62 mg of <b>AE-382-2</b> was transferred into the vial of the 5 mL Schlenk flask in air. Afterwards the flask was set under argon.	<b>AE-382-MS</b>

	14:00	<p>The sample was submitted to Nicole Fritz. Further work done by Nicole:</p> <p>Rough description of sample handling: The flask was opened to air and a small amount (approx. 0.1 mg) was taken out and diluted in acetonitril (ACN) 3 times. After taking out the small amount the flask was closed again. Afterwards the sample was injected into the MS device. Before and after the measurment a calibration measurement was performed. (<b>AE-382-MS ACN_1</b>)</p> <p>After approx. 10 min a sample of the same solution was measured (<b>AE-382-MS ACN_1_10min</b>)</p> <p>A new sample was prepared the same way as the first sample (<b>AE-382-MS ACN_2</b>)</p> <p>A new sample was prepared analog to the first sample, but the small amount taken out was first dissolved in water (approx. 0.5 mL) and than diluted with ACN (<b>AE-382-MS H2O_ACN</b>)</p>	
31.03	/	The elemental analysis of <b>AE-382-EA</b> was measured by Sandra Köhn	

## Analysis

Date	Time	Sample name	Analysis method	Used device	Solvent	Raw Data	Processed Data	Used for plotting	Comparative Data	Interpretation
18.11	13:21	AE-382-1	NMR 1H, 31P{1H}quant	IAAC 400 MHz I	DMSO-d6	AE-382-1.zip	AE-382-1_10.nmrium	/	AE-310: Synthesis of [Ru(CO)(OH)2(PNPtBu)] using Ag2O, 5 h	As expected. Only minor impurities
12.12	/	AE-382-3	NMR 1H, 31P{1H}quant	IAAC 400 MHz I	DMSO-d6	AE-382-3.zip	2024-12-13_12-58-51.zip_AE-382-3_10.nmrium	/	/	As expected. Only minor impurities
15.01	17:31	AE-382-2	NMR 1H, 31P{1H}quant	IAAC 400 MHz I	MeCN-d3	AE-382-2.zip	AE-382-2_10.nmrium	/	/	As expected. Only minor impurities
13.03	14:10	AE-382-MS ACN_1	ESI-MS	ESI-TOF-Bruker	MeCN	AE-382-MS-ACN_1.xy	AE-382-2-MS-ACN_1.txt AE-382-MS-ACN_1.png	ESI-MS.py	C24H43NO3P2RuCO.csv calculated via <a href="https://www.envipat.eawag.ch/index.php">https://www.envipat.eawag.ch/index.php</a>	Very weak singal which could be M+H, but significant signals at higher mass, full interpretation pending Main signal could be from carbonate forming and removing both OH groups through addition of CO2
13.03	14:20	AE-382-MS ACN_1_10min	ESI-MS	ESI-TOF-Bruker	MeCN	AE-382-MS-ACN_1_10min.xy	AE-382-2-MS-ACN_1_10min.txt AE-382-MS-ACN_1_10min.png	ESI-MS.py	C24H43NO3P2RuCO.csv calculated via <a href="https://www.envipat.eawag.ch/index.php">https://www.envipat.eawag.ch/index.php</a>	Significantly different pattern -> most likely decomposition
13.03	14:30	AE-382-MS ACN_2	ESI-MS	ESI-TOF-Bruker	MeCN	AE-382-MS-ACN_2.xy	AE-382-2-MS-H2O_ACN.txt AE-382-MS-ACN_2.png	ESI-MS.py	C24H43NO3P2RuCO.csv calculated via <a href="https://www.envipat.eawag.ch/index.php">https://www.envipat.eawag.ch/index.php</a>	Very weak singal which could be M+H, but significant signals at higher mass, full interpretation pending Weaker signal than ACN_1. signal at 627 less pronounced. Main signal at 586 u (delta to expected: 28 u)
13.03	14:35	AE-382-MS H2O_ACN	ESI-MS	ESI-TOF-Bruker	MeCN + H2O	AE-382-MS-H2O_ACN.xy	AE-382-2-MS-ACN_2.txt AE-382-MS-H2O_ACN.png	ESI-MS.py	C24H43NO3P2RuCO.csv calculated via <a href="https://www.envipat.eawag.ch/index.php">https://www.envipat.eawag.ch/index.php</a>	Very weak singal which could be M+H, but significant signals at higher mass, full interpretation pending Some new, high mass peaks
31.03	/	AE-382-EA	elemental analysis	?	/	AE-382-EA.jpg	C: 49.78 % H: 6.94 % N: 2.45 %	/	C: 51.6 % H: 8.12 % N: 2.51 %	C and H contant to low by 1.8 pp for C and 1.2 pp for H. N content fits roughly

## Product characterization and results

Sample	Mass [mg]	Purity	Mass <sub>pure</sub> [mg]	Amount [μmol]	Yield [%]	Description	Storage location
AE-382-1	149.79	> 96 %	/	0.268	32	grayish olive-ochre very fine, sticky solid, <a href="#">Product.jpg</a>	used
AE-382-3	40.87	> 96 %	/	/	/	grayish olive-ochre very fine, sticky solid	freezer drawer A
AE-382-2	23.87	> 96 %	/	/	/	grayish olive-ochre very fine, sticky solid	freezer drawer A

MS fits to +CO -HO.

## Linked experiments

- [AEI-010: Synthesis of \[Ru\(CO\)\(OH\)2\(PNN\)\] using Ag2O](#)

- [JSC-KS-04: Degassing of n-Heptane](#)

- [AEI-024: Synthesis of \[Ru\(CO\)\(OH\)2\(NNN\)\] with 2,6-bis\(N,N-diethylaminomethyl\) using Ag2O](#)

- [AEI-026: Capillary production](#)

- [AEI-048: Synthesis of \[Ru\(CO\)\(OH\)2\(NNN\)\] with 2,6-bis\(N,N-diisopropylaminomethyl\) using Ag2O](#)

- [AEI-064: Synthesis of \[Ru\(CO\)\(OH\)2\(PNP\)\] with 2,6-bis\(P,P-di\(tertbutyl\)phosphinomethyl\) using Ag2O](#)

- [AEI-073: Synthesis of \[Ru\(CO\)\(OH\)2\(PNP\)\] with 2,6-bis\(P,P-di\(tertbutyl\)phosphinomethyl\) using Ag2O](#)

- [AE-101: Synthesis of \[RuCl2\(CO\)\(PNP\)\] with 2,6-Bis\(\(tert-butyl\)phosphinomethyl\)pyridine via \[RuCl2\(CO\)\(p-cymene\)\]](#)

- [AE-102: Degassing of THF](#)

- [AE-105: Synthesis of \[Ru\(CO\)\(OH\)2\(PNP\)\] with 2,6-bis\(P,P-di\(tertbutyl\)phosphinomethyl\) using Ag2O](#)

- [AE-108: Synthesis of \[RuCl2\(CO\)\(PNP\)\] with 2,6-Bis\(\(tert-butyl\)phosphinomethyl\)pyridine via \[RuCl2\(CO\)\(p-cymene\)\]](#)

- AE-114: Synthesis of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$  with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using  $\text{Ag}_2\text{O}$
- AE-121: Degassing of n-Heptane
- AE-145: Synthesis of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$  with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using  $\text{Ag}_2\text{O}$
- AE-152: Degassing of MeCN
- AE-157: Synthesis of  $[\text{RuCl}_2(\text{CO})(\text{PNP})]$  with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via  $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$
- AE-167: Irradiation of PhPDA (365 nm + white LED, gas phase measurement), AE-151, 1.5 mg/mg SDS, 2 mg/mL PhPDA
- AE-198: Synthesis of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP}\{\text{tBu}\})]$  with  $\text{Ag}_2\text{O}$  with different reaction times (1, 2, 5, 24 h)
- AE-KRA-204: Synthesis of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$  using  $\text{Ag}_2\text{O}$ , 5 h
- AE-208: Degassing of Milli-Q water
- KRA-080: Degassing of THF
- KRA-104: Degassing of Milli-Q water
- AE-230: Synthesis of  $[\text{RuCl}_2(\text{CO})(\text{PNP})]$  with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via  $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$
- AE-232: Synthesis of  $[\text{RuCl}_2(\text{CO})(\text{PNP})]$  with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via  $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$
- AE-238: Synthesis of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$  using  $\text{Ag}_2\text{O}$ , 5 h
- AE-246: Radiation of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP}\{\text{tBu}\})]$  (AE-238, 97 % pure) 365 nm and white
- AE-252: Synthesis of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$  using  $\text{Ag}_2\text{O}$ , 5 h
- AE-261: Synthesis of  $[\text{RuCl}_2(\text{CO})(\text{PNP})]$  with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via  $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$
- AE-281: Stability test of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP}\{\text{tBu}\})]$  (AE-252) under air and afterwards irradiation with 365 nm and white
- AE-310: Synthesis of  $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$  using  $\text{Ag}_2\text{O}$ , 5 h



## Linked resources

Equipment - [Freezer Lab E004 CEEC II](#)

Protocol - [Preparation of NMR Sample](#)

Protocol - [Filtration with frit technique](#)

Protocol - [Schlenk Technique](#)

Protocol - [Preparation of capillary for H<sub>2</sub>O-NMR](#)

## Attached files

AE-382-EA.jpg

sha256: 4f383ed83f0808ff6e6830a4ea048a0522ef4c7b5a5ec1265b0a43b53c06d46f



C24H43NO3P2RuCO.csv

sha256: e324012daf18017a2ef1d1d4b926e79458494cebb8e25ae6d29361227fe719fd

ESI-MS.py

sha256: c5a21c33581840e3dd1f8fc8697c254996c7dc42ae6121e877fb052a499b9e58

AE-382-2-MS-H<sub>2</sub>O\_ACN.txt


sha256: abf859520850a5b424900e4f39f7c0f01901a7becbe0999eb2ca0f9091f33221

AE-382-2-MS-ACN\_2.txt

sha256: ab0cfb93750befe6ae85bff2c668ed879a7a4319ff120611cdbe4db5275721bb

AE-382-MS-H<sub>2</sub>O\_ACN.png

sha256: f40efc75601edc7db4941c312b70af01d2188cad98a60f521b4b24477eb6243



AE-382-2-MS-ACN\_1.txt

sha256: 424db80972ca417377214c9cb93c1a6b5b91659ae3d9a9ff20d1d670d9fe20f8

AE-382-MS-ACN\_1\_10min.xy


sha256: 676c891e91704f16c7612aaf3389ac1b7207f38de6e7a812e0094d828aac86dd

AE-382-MS-ACN\_1.xy

sha256: 02820f8babd2e18759aaef8daa061ac389943db264f20d88e7a483bb3161b2bb


AE-382-MS-ACN\_1.png

sha256: 5c0938c4b182ca1c8e457c827fa725597981118f836ee69e5c4e1a872c1daf77



AE-382-MS-ACN\_1\_10min.png

sha256: ea9d1defa57854a2c424ab801c93f74fc594393b8dd174c5b16fa56a3e136a16




AE-382-2-MS-ACN\_1\_10min.txt

sha256: f1b847c4acba05e3bd3e688bd266c0f888576d67c2c09cc6d72d6dca1590330d

AE-382-MS-ACN\_2.png

sha256: ceb94ce6feb1e1753ee7495f08c2c6e2a5a570b28ce6b46cab171a06dcb2d864



AE-382-MS-ACN\_2.xy


sha256: f9d28caeff5eb4cb00119f745d2dfe0cf7f46deb86d4206bc05c16ff8858b750

AE-382-MS-H2O\_ACN.xy

sha256: 856d9846489a7abe2b0de70ce7464150c2ae07fdced95c38c02a11a9a9508253

AE-382-2.jpg

sha256: 38f8235762d0898057456031830bf15ea796c9815115412088ace58e177c60e8



AE-382-2\_10.nmrium

sha256: dd67522ebeda93741474bc5723fed107f200e7270475dce713868e9864864ca6

AE-382-2.zip

sha256: 0a796ce852827145d822c2801f3a9a793c276888db3c0b09c96f0abcc4e39559

AE-382-3.zip


sha256: 0cd0258b2d0b5acd6e1d5d585a472b5b21b1bc4e44b502a3764c0e607cd43c79

2024-12-13\_12-58-51.zip\_AE-382-3\_10.nmrium

sha256: 2c1adc0bc9ecfa25d8678ebc1e96697b01aff93b56d17b7be86af41934ee0772

3-in-DMSO-d6.jpg

sha256: dd6952fbdcc9b2b8e6c171ce96e999d58c6b850d69866d35a669816d0ef32347



AE-382-3.jpg

sha256: af7c23252fae8d6e1e9056e27eb399d4fbb91ecaad6b624c7332a5a8f3685880



before-filt.jpg

sha256: 0bd435be27f7cb56ea398bb8f2c39c2bd65134acec43efe7b55a0bb97b14710f




in-THF.jpg

sha256: b10d22eb0f942989ac32e1379099c1693177a02a2dd919296b38e135be908a86



after-n-hept.jpg

sha256: c206c51b6df0041594e2dc1483adc2ef125658cb21b9e4714e1f87e7e125d400



washing-of-filtrate-on-1st-frit.jpg

sha256: f8bf5c7f8ea54f2928ed91255bbf3b1457519138ee013ccda3cac7218fe06dd8



AE-353.png

sha256: cfc2963bdd5bfa9c84d644f2a2268bb1298a88bd362365e11a2098efe3173a77



AE-382-1\_10.nmrium

sha256: 105e1c8e653ddd70ce6e280e15fc0654b8480fc2a1e0c5d596bb9c24f0865d41

AE-382-1.zip

sha256: ae94ad6033d2a53eb8fef04b849b9fac631f5bd96b749669626d46b132eae5e

after-reaction.jpg

sha256: f16e4755b79d5a5939b25cc5e1761b3852cacfbfb3609e6120ba499015ae358d



Product.jpg

sha256: 6e4f492beba5706fb904ece1f5979d37a1fac8513434f0091a68bdbba2857053



after-adding-n-hept.jpg

sha256: 9081de487db49b9b43a8901601b95078f05b69c3478efecd764d4316755e8aa4




Start-of-reaction.jpg

sha256: e4e2a323423e85b1b2e68b98b457541f201890bd7179c807063d99e61c3b2916



AE-382-1.jpg

sha256: 555cb7dce37d606c6ff3860d0c15a178bf5a4872e5d66b505f55f82813c9046e



after-adding-THF.jpg

sha256: 2b3b591b6e78597aac125ed4fb6cb1ca996867bd973e1e4fbb7831c338bea000



Solid-after-night.jpg

sha256: 80072f5c2fb15e55a786a9b31a9883328501a0e964a7fc1639d3ae5a6beb43e4



AE-20421.cdxml

sha256: 70d2a67f6ae7477a3c31ec8e776e001181c617b18e0d276cd2995bdbdfdb2ca2

## Comment

On 2024-12-13 13:41:11 Jacob Schneidewind wrote:

Notes:

\* adding future recommendation to change MeCN purification step (same as AE-302-1)



Unique eLabID: 20241113-5686e080afe9a3b141d8cfff3b0f1b849e0097d

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1502>