

JU-KRA-172: Degassing of Diethylether

Date: 2024-08-19
Tags: Degassed solvent KRA JU
Status: Done
Created by: Jens Ulbrich

Reagents

Name	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
Diethylether (SPS MS Braun Lab 107 ZAF)	/	/	/	/	/	approx. 300 mL
Molecular sieve (Experiment - AE-291: Activation of molecular sieve)	/	/	/	/	/	/

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard Schlenk technique or in an argon filled glovebox (MBRAUN).

Date	Time	Step	Observations
19.08.24	10:30	J.Young schlenk flask was prepared	
		Molar sieve was added to the flask (approx. 1/4)	* took to long due to the bigger size of the molecular sieve
		Collecting at SPS according to Protocol - Collecting of SPS solvents	
	15:29 - 15:49	The solvent was degassed according to Protocol - Degassing of organic solvents	

Linked experiment

- [AE-291: Activation of molecular sieve](#)

Linked resources

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [KFT \(Karl-Fischer-Titration\)](#)

Protocol - [Collecting of SPS solvents](#)



Unique eLabID: 20240819-270969cf6034b16f10d838f29f02def88e6a861a

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1293>