

AEI-046: Degassing of NEt₃ (Triethylamin)

Date: 2022-11-29
Tags: AEI | Degassed
Status: Done
Created by: Alexander Eith

Reaction scheme/sample structure

Reagents

Name	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
Triethylamin						approx. 15

Procedure/observations

Date	Time	Step	Observations
29.11	13:20	Triethylamin was degassed according to [Protocol] Freeze Pump Thaw . The liquid was stored over activated molar sieve (3A)	

Analysis

Date	Time	Sample name	Analysis method	Solvent	File name	Interpretation

Product characterization

Linked resources

Protocol - [Freeze Pump Thaw](#)

Protocol - [Degassing of fluids in normal Schlenk flask](#)



Unique eLabID: 20221129-b1536ea88c8d5d491e76e718b96ed11386d247ca

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=289>