

JSC-KS-20: Degassing of MeOH

Date: 2022-12-02

Tags: ☐ Degassed ☐ JSC-KS

Status: Done

Created by: Katrin Selzner

Reaction scheme/sample structure

Reagents

Name	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
Methanol						approx. 500

Procedure/observations

Date	Time	Step	Observations
01.12.22		A 500 mL-solvent flask was prepared. (AEI)	
		Molar sieve 3 A (approx. 100 mL) was filled in the flask and the atmosphere was subsequently changed to argon.	
01.12	14:00 - 14:25	THF was transferred into the flask according to [Protocol] Cannula Filtration/Transfer	
	14:25 - 15:00	Argon was sparged through the solution analog to [Protocol] Degassing of organic solvents for 30 min	
06.12	9:45	The water content was determined according to [Protocol] KFT (Karl-Fischer-Titration)	322.7 mg sample
22.12	9:00	KFt according to [Protocol] KFT (Karl-Fischer-Titration)	150.8 mg

Analysis

Date	Time	Sample name	Analysis method	Solvent	File name	Interpretation
06.12	9:45		KFT			Amount H2O: 6.8 µg, 21.1 ppm
22.12	9:00		KFT			Amount H2O: 3.3 µg, 21.9 ppm

Product characterization

Linked resources

Protocol - [Cannula Filtration/Transfer](#)

Protocol - [Degassing of fluids in normal Schlenk flask](#)

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [KFT \(Karl-Fischer-Titration\)](#)



Unique eLabID: 20221202-8c116618382d22d5b6503d549260b5e6f3072712
Locked by Katrin Selzner on 2023-02-06 at 17:38:18
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=306>