

AE-396: Degassing of THF

Date: 2024-12-05

Tags: Degassed AE

Category: Prep work

Status: Done

Created by: Alexander Eith

Reagents

Name	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
THF	/	/	/	/	/	approx. 450 mL

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard Schlenk technique or in an argon filled glovebox (MBRAUN).

Date	Time	Step	Observations
04.12	14:00	THF was collected from the SPS according to Protocol - Collecting of SPS solvents	
	14:35 - 15:08	The THF was transferred into a Young Type-Schlenk flask and degassed according to Protocol - Degassing of organic solvents starting from step 5	Degassing was done for 30 min

Analysis

Date	Time	Sample name	Analysis method	Solvent	Raw Data	Processed Data	Interpretation
------	------	-------------	-----------------	---------	----------	----------------	----------------

Linked experiment

- [AE-289: Degassing of THF](#)

Linked resources

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [KFT \(Karl-Fischer-Titration\)](#)

Protocol - [Collecting of SPS solvents](#)



Unique eLabID: 20241205-23a7dbd7b5501912681b4a471276ed1a2fe1bbca

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1555>