

AE-442: Degassing of Milli-Q water

Date: 2025-02-06
Tags: Degassed AE
Category: Prep work
Status: Done
Created by: Alexander Eith

Reaction scheme/sample structure

Reagents

Name	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
Milli-Q water	/	/	/	/	/	approx. 90 mL

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard Schlenk technique or in an argon filled glovebox (MBRAUN).

Date	Time	Step	Observations
05.02		A 100 mL young type flask was prepared	
		Milli-Q water was added	
	11:20 - 11:51	Argon was sparged through the solution analog to [Protocol] Degassing of fluids for 30 min	

Analysis

Date	Time	Sample name	Analysis method	Solvent	Raw Data	Processed Data	Interpretation

Product characterization

Linked experiment

- [AE-208: Degassing of Milli-Q water](#)

Linked resources

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [KFT \(Karl-Fischer-Titration\)](#)

Protocol - [Collecting of SPS solvents](#)



Unique eLabID: 20250206-13ed8caa91c4467777ae8d807cd1b365baddbd88
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1712>