

AE-460: Synthesis of ^{17}O labeled $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$

Date: 2025-03-12


Tags: Dihydroxo Ag 20 PNP
P(tBu)N(py)P(tBu) $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$
NMR AE Future reference analytics ^1H
 ^{31}P ^{17}O

Category: Two-photon

Status: Done


Created by: Alexander Eith

Reaction scheme/sample structure



Chemical Formula: $\text{C}_{24}\text{H}_{43}\text{Cl}_2\text{NOP}_2\text{Ru}$
Molecular Weight: 595,53152

$[\text{RuCl}_2(\text{CO})(\text{PNP})]$



Chemical Formula: $\text{C}_{24}\text{H}_{45}\text{NO}_3\text{P}_2\text{Ru}$
Molecular Weight: 558,64552

$[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP})]$

Silver Oxide

Molecular Weight:
231.74 g/mol

Literature/reference experiments

Literature	https://doi.org/10.1039/D1EE01053K
Reproduction	/
Related experiments	Two-photon - AE-382: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag 20 , 5 h Two-photon - AE-456: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag 20 , 5 h

Reagents

Name	CAS-Number/Experiment-Number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume _{theo} [ml]	Volume _{exp} [ml]	Density [g/mL]
$[\text{RuCl}_2(\text{CO})(\text{PNP}\{\text{tBu}\})]$	Two-photon - KRA-DS-02: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(p\text{-cymene})]$ I	0.118	1.00	70.47	70.76	595.53	/	/	/
Ag ^{17}O (10 % ^{17}O)	Two-photon - AE-450: Synthesis of Ag ^{17}O from AgNO $_3$ and H $_2$ (^{17}O)	0.727	6.14	168.38	167.60	231.74	/	/	/

H ₂ ¹⁷ O (10 % 17O)	Two-photon - AE-459: Degassing of H ₂ (17)O	27.69	234	499	/	18.02	0.500	0.50	0.998
THF	Prep work - AE-396: Degassing of THF	/	/	/	/	/	3.5	3.6	/

Work-up Reagents

Name	CAS Number / Experiment-Number	Inventory number	Mass _{exp} [g]	Volume [ml]	Concentration [M]
THF	Prep work - AE-396: Degassing of THF	/	/	7.5	/
n-heptane	AE-315: Degassing of n-Heptane	/	/	14 + 3	/
MeCN-d ₃	KRA-056: Degassing of MeCN-d ₃	/	/	1	/

Procedure/observations

All steps, unless mentioned otherwise, were carried out under argon using standard Schlenk technique or in an argon filled glovebox. (see Protocol [Protocol - Schlenk Technique](#))

Due to the instability of Ag₂O and Dihydroxo to the light, the procedure was done mostly in the dark, by either wrapping the flasks containing either or both in Aluminium foil or by tuning off the light in the lab and fumehood

Date	Time	Step	Fotos	Observations
12.03		A 25 mL brown glas Schlenk flask was prepared, by heating under vacuum	/	/
	9:10	[RuCl ₂ (CO)(PNP)] and Ag ₂ (17)O was weighted in air into the reaction schlenk flask.	/	/
		The light in the fumehood was turned out	/	/
	9:25	The flask was set under argon	/	/
	9:42	THF (3.6 mL) was added to the reaction schlenk flask.	/	/
	9:45	H ₂ (17)O (0.50 mL) was added to the reaction schlenk flask	/	/
	9:45 - 14:45	The reaction mixture was stirred for 5 h at 500 rpm. During this the flask was wrapped with aluminium foil	/	/
		In the meantime a 100 mL Schlenk flask was prepared	/	/
	14:45 - 14:55	The reaction mixture was allowed to settle	/	/
		<i>Work-up:</i>		
	14:55	The obtained suspension was filtered through a PTFE syringe filter (pore diameter: 0.2 µm) in the new 100 mL flask	after filt through PTFE.jpg	purplish solution

	15:30 - 15:32	The solvents of the filtrate were removed under reduced pressure, for this no water bath was used	/	The mixture solidified
	15:33	A water bath (rt) was put under the flask, while keeping the flask under an reduced pressure	/	/
	15:40 - 16:18	The obtained solid was dried under reduced pressure	after drying.jpg	AE-460-0 , brown film
		The liquid in the cooling trapped was transferred into a 10 mL snap on cap vial. It was stored under air in the SSC at rt	liquid from drying.jpg (mislabeld as -2 in Photo)	AE-460-1
	16:20	AE-460-0 was stored in the dark at -20 °C under argon	/	/
		From now on AE-460-0	next morning.jpg	no change
13.03	10:10	To the solid THF (7.5 mL) was added und the obtained mixture was stirred at 500 rpm. In between the flask was shaken for approx. 30 s.	in THF.jpg	light brwon, pinkish suspension
	10:17	n-heptane (14 mL) was added to the mixture under stirring	addition of heptane.jpg setteled.jpg	pinkish solution, yellowish solid
	10:27 - 10:54	The mixture was filtered through a pore 4 frit according to Protocol - Filtration with frit technique into a 100 mL-Schlenk	during filt.jpg	no change
	10:54 - 11:00	The solid was washed with n-heptane (3 * 1.0 mL)	product.jpg	brwonish solid
	11:02 - 13:05	The solid was dried under reduced pressure	/	AE-460-2
		2 10 mL Schlenk flask (1 with a stir bar, 1 with a 2 mL screw cap vial) were prepared	/	/
		From now on AE-460-2:		
	13:30 - 13:40	The solid was transferred into a 2 mL vial in a 10 mL Schlenk flask in air, thereby it was weighed	/	/
		During this an aliquote (8.97 mg) was transferred into the other 10 mL Schlenk flask	/	/
		Both flask were set under argon	/	/
		AE-460-2 was stored in Equipment - Freezer Lab E004 CEEC II in the dark at -20 °C under Ar	/	/
	15:28	To the aliquote MeCN-d3 (1 mL was added)	solution in MeCN-d3.jpg	orange brown solution
	- 15:32	The obtained solution was stirred at 600	/	/
	15:33	From this an NMR sample was prepared according to Protocol - Preparation of NMR Sample during this the solution was filtered through an PTFE syringe filter (pore diameter: 0.2 µm)	/	/

	15:45	The NMR was submitttd for an 1H NMR	NMR sample.jpg	AE-460-2-1 yellow orange solution
14.03	12:00	The same NMR tube was submitted for an 17O NMR	/	AE-460-2-2
19.03	11:10	The same NMR tube was submitted for several NMR experiments		AE-460-2-3
26.03	16:00	The same NMR tube was submitted for 13C NMR experiment on 500 MHz device	AE-460-2-4.jpg	AE-460-2-4 very small amount of solid observed at the wall, otherwise no change, by eyesight should be still ok

Analysis

Date	Time	Sample name	Analysis method	Used device	Solvent	Raw Data	Processed Data	Comparative Data	Interpretation
13.03	16:01	AE-460-2-1	NMR 1H	ZAF 300 MHz OA	MeCN-d3	AE-460-2-1.zip	AE-460-2-1_10.nmrium	Two-photon - AE-382: Synthesis of [Ru(CO)(OH)2(PNpTbU)] using Ag2O, 5 h (AE-382-2)	As expected. Very few impurities observed.
14.03	12:46	AE-460-2-2	NMR 1H, 31P{1H}, 17O{1H}	IAAC 400 MHz	MeCN-d3	AE-460-2-2.zip	AE-460-2-2_10.nmrium	Two-photon - AE-382: Synthesis of [Ru(CO)(OH)2(PNpTbU)] using Ag2O, 5 h (AE-382-2)	In 1H and 31P as expected. Very few impurities observed. In 17O: 3 signals observed. As expected all broad singulets. Shifts (intensity): 19.8 (43 %), 94.7 (47 %) and 244.0 (9.5 %) ppm. Main signals would roughly fit to OH groups (expected shift range for alcohols: very approx. -50 to 100 ppm according to https://chem.ch.huji.ac.il/nmr/techniques/1d/row2/o.html) Signal at 244 ppm: no good idea, shift could fit to peroxide, carbonate or carboxylic acid according to https://chem.ch.huji.ac.il/nmr/techniques/1d/row2/o.html
19.03	15:36	AE-460-2-3	NMR 1H, 13C{1H}, 31P{1H}, 13C{1H}-APT, 1H/1H-COSY, 1H/13C{1H}-HMBC, 1H/13C{1H}-HSQC, 1H/31P{1H}-HMBC	IAAC 400 MHz I	MeCN-d3	AE-460-2-3-all-spectra.zip AE-460-2-3-1H.zip AE-460-2-3-31P.zip AE-460-2-3-1H/31P-HMBC.zip AE-460-2-3-13C.zip AE-460-2-3-1H/13C-HSQC.zip AE-460-2-3-1H/13C-HMBC.zip AE-460-2-3-1H-COSY.zip AE-460-2-3-13C-APT.zip	AE-460-2-3-1H-13C-31P.nmrium AE-460-2-3-1H-13C-31P-comparison.nmrium AE-460-2-3-31P-HMBC.nmrium AE-460-2-3-13C-HSQC.nmrium AE-460-2-3-13C-HMBC.nmrium AE-460-2-3-1H-COSY.nmrium	AE-460-2-2	1H, 31P as expected, all signals identified 13C: Evaluation pending, some triplets (?) observed around 38 ppm 31P-HMBC: all cross signals observed 13C-HSQC: all cross singals observed 13C-HMBC: evaluaion pending, broad couplings 1H-COSY: Cross signals observed for aromatic signals and for methylene signals
28.03	12:56	AE-460-2-4	NMR 1H, 13C	IOMC 500 MHz	MeCN-d3	AE-460-2-4.zip	AE-460-2-4_101.nmrium	AE-460-2-2 AE-460-2-3	13C: all signals observed, signals of CH2 and Cq of tBu no triplett observed 1H: mainly as expected, but for tBu groups more complex coupling patterns observed. Looks like weak douplet of tripletts

Product characterization and results

Sample	Mass [mg]	Purity	Mass _{pure} [mg]	Amount [μmol]	Yield [%]	Description	Storage location
AE-460-2	26.51	>98 %	/	47.5	40	dried product.jpg light greenish gray solid	Equipment - Freezer Lab E004 CEEC II
Results	High yield obtained. Singlas observed in 17O, two of the 3 would fit to dihydroxo groups						

Future recommendations

Old procedure	Problem	Suggested new procedure
---------------	---------	-------------------------

/	/	always use self made silver(I) oxide
---	---	--------------------------------------

Linked experiments

- AEI-010: Synthesis of [Ru(CO)(OH)₂(PNN)] using Ag₂O
- JSC-KS-04: Degassing of n-Heptane
- AEI-024: Synthesis of [Ru(CO)(OH)₂(NNN)] with 2,6-bis(N,N-diethylaminomethyl) using Ag₂O
- AEI-026: Capillary production
- AEI-048: Synthesis of [Ru(CO)(OH)₂(NNN)] with 2,6-bis(N,N-diisopropylaminomethyl) using Ag₂O
- AEI-064: Synthesis of [Ru(CO)(OH)₂(PNP)] with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag₂O
- AEI-073: Synthesis of [Ru(CO)(OH)₂(PNP)] with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag₂O
- AE-101: Synthesis of [RuCl₂(CO)(PNP)] with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via [RuCl₂(CO)(p-cymene)]
- AE-102: Degassing of THF
- AE-105: Synthesis of [Ru(CO)(OH)₂(PNP)] with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag₂O
- AE-108: Synthesis of [RuCl₂(CO)(PNP)] with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via [RuCl₂(CO)(p-cymene)]
- AE-114: Synthesis of [Ru(CO)(OH)₂(PNP)] with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag₂O
- AE-121: Degassing of n-Heptane
- AE-145: Synthesis of [Ru(CO)(OH)₂(PNP)] with 2,6-bis(P,P-di(tertbutyl)phosphinomethyl) using Ag₂O
- AE-152: Degassing of MeCN
- AE-157: Synthesis of [RuCl₂(CO)(PNP)] with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via [RuCl₂(CO)(p-cymene)]
- AE-167: Irradiation of PhPDA (365 nm + white LED, gas phase measurement), AE-151, 1.5 mg/mg SDS, 2 mg/mL PhPDA
- KRA-056: Degassing of MeCN-d₃
- AE-198: Synthesis of [Ru(CO)(OH)₂(PNP{tBu})] with Ag₂O with different reaction times (1, 2, 5, 24 h)

- AE-KRA-204: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h

- AE-208: Degassing of Milli-Q water

- KRA-080: Degassing of THF

- KRA-104: Degassing of Milli-Q water

- AE-230: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$

- AE-232: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$

- AE-238: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h

- AE-246: Radiation of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP}\{\text{tBu}\})]$ (AE-238, 97 % pure) 365 nm and white

- AE-252: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h

- AE-261: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$

- AE-281: Stability test of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNP}\{\text{tBu}\})]$ (AE-252) under air and afterwards irradiation with 365 nm and white

- AE-310: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h

- AE-315: Degassing of n-Heptane

Prep work - AE-396: Degassing of THF

Prep work - AE-442: Degassing of Milli-Q water

Two-photon - AE-382: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h

Two-photon - AE-450: Synthesis of $\text{Ag}_2(17)\text{O}$ from AgNO_3 and $\text{H}_2(17)\text{O}$

Two-photon - KRA-DS-02: Synthesis of $[\text{RuCl}_2(\text{CO})(\text{PNP})]$ with 2,6-Bis((tert-butyl)phosphinomethyl)pyridine via $[\text{RuCl}_2(\text{CO})(\text{p-cymene})]$ I

Two-photon - AE-456: Synthesis of $[\text{Ru}(\text{CO})(\text{OH})_2(\text{PNPtBu})]$ using Ag_2O , 5 h

Two-photon - AE-459: Degassing of $\text{H}_2(17)\text{O}$

Linked resources

Equipment - [Freezer Lab E004 CEEC II](#)

Protocol - [Preparation of NMR Sample](#)

Protocol - [Filtration with frit technique](#)

Protocol - [Schlenk Technique](#)

Protocol - [Preparation of capillary for H₂O-NMR](#)

Attached files

AE-460-2-4_101.nmrium

sha256: 0f0f0decb6f8f4057faba0975df3538e51e991dc415ceb6ddcc6f416e0c371b4

AE-460-2-4.zip

sha256: 5b7fbb0bd0dd7f23c1fe8ec0ca126f1871eb0dac74996f01142c7e1245b33c0c

AE-460-2-4.jpg

sha256: 339cf80804eeb31086c46d2f768852271b67b2a73780edf501f4a2dd3773c1ca



AE-460-2-3-1H-COSY.nmrium

sha256: b6a13a70253a959237109b9b3c133f660abf7d01032468da12a1a8c6b95d0e3f

AE-460-2-3-13C-HMBC.nmrium

sha256: 0693dd6337b7b1033285a4b677856db0cc6370c8794e5f741afb3496f991ba1c

AE-460-2-3-13C-HSQC.nmrium

sha256: f15c877495636d780c6850a74be42fb15360931c4f7ab39369ee8133ce8ba622

AE-460-2-3-31P-HMBC.nmrium

sha256: e863b8b93ee216de577000b6aae787d443f6198fe59feee59822e5205bd30933

AE-460-2-3-1H-13C-31P-comparison.nmrium

sha256: 387723c312249207b374cca84dc499b0701233c16f917579b1602635211e5111

AE-460-2-3-1H-13C-31P.nmrium

sha256: 2861705653fabb933bf1c658e3dacad5e56a22cc1199359ccc073645275373da

AE-460-2-3-13C-APT.zip

sha256: 1346a8f6e9aea0d724ed317d8f571445b641269428219dbdfa79d5fad2a38a5b

AE-460-2-3-1H1H-COSY.zip

sha256: 755604ec50ebf08c12661834ceb4676dd9fcd3e50bf30d6edb67f4f31d8f826c

AE-460-2-3-1H13C-HMBC.zip

sha256: 2cac99b395ca72170ced55e76c716c2e81ddd2440051cafd3e6c064452ae9c90

AE-460-2-3-1H13C-HSQC.zip

sha256: 563bedf76890d716040aa6ed94ac26fd9a4c90e6889b444574ffba68007d77e1

AE-460-2-3-13C.zip

sha256: c91aa0b710f09214bba6d96daf34a744db69bcd460fb3e9ba64ecf040947ebf

AE-460-2-3-1H31P-HMBC.zip

sha256: 931a9517b9db70f6fcd38e22a8b30d258e21037923d856a6d9bd5a540d64d2cf

AE-460-2-3-31P.zip

sha256: fc62cd7077fcb1b55b791f53212433de096087ce999dbf314284ceb180baa139

AE-460-2-3-1H.zip

sha256: e535ca68a917abffba9f224a48a0592bf121ba262ba2b4f66896279887c03b89

AE-460-2-3-all-spectra.zip

AE-460-2-2_10.nmrium

sha256: 38d9b61bcb291328304cc823db02a61e34cdc5d2afce8bba5282c6ae677611c2

AE-460-2-2.zip

sha256: 7e5d9fbeeec1351960301847f5b39e4a6d6f545d94ec3ca4d31b94c3e33e418b7

NMR-sample.jpg

sha256: 550d45692d2125ef63317317e9ff880f51ba7731d1740617c397b9585711b023



solution-in-MeCN-d3.jpg

sha256: 81258d0996ed6c866b974705e3f0e97d356679fb166aea0f045bd0af76253201



dried-product.jpg

sha256: 5ab558cf962bbfab351bcdbedb0d90c6cff1eacfd47e14d6bd8a47dee34e0eb



product.jpg

sha256: db8f20b9be24b7b22b692500850f24ccc30463dcf835702992201593751e9275



during-filt.jpg

sha256: d67b7cf6eb21fa2764b91d08699681ba03bacb57a94aa8066962e15c5a924a20



addition-of-heptane.jpg

sha256: e28b7fa7c89dba0b67f577f80544e5f28dfaf5c4ac8e0cb5d2f53144e8ec657a



in-THF.jpg

sha256: 6d56e3f739d12e5d8228ad551d43b3fa237aee63391827938495a626c51a8ce4



setteled.jpg

sha256: b2c6dfcf2a9a2a109a0ba197af562cb6999ac90121fbb31c64feba7f4fb859e2



next-morning.jpg

sha256: 96585dbfe4c65b1fd65daa9439bf894c638c4d25e736a8e139d3da69ee78b945



liquid-from-drying.jpg

sha256: 4dfe83cf6e193ccc094be42dd863e4cec22aab178080acf76b57bee793a8a611




after-filt-through-PTFE.jpg

sha256: b7da4db13440dcdffe339519cbc73c3c9d602789af884f59b68805a3bb7e4bff



after-drying.jpg

sha256: 19a524a9130687779a43e3732bcaf5e6d4a121d7b2808c9db07561cf27e14770



AE-460-2-1_10.nmrium

sha256: dc41cd158c181a0233bd57f5214d0ce2e75d9cad7e2ee423af404e4a2078b351

AE-460-2-1.zip

sha256: c6c349925d0fc577de3fd2fbff460a65b99e7300c3a8804d3b0ab2edd4903c2b

AE-460.png

sha256: a631377ee73b6a2b721398cf2e357b33404e7647968a7fbdbb6035fc536dedd2



AE-460.cdxml

sha256: 7081e3701fedd9c03b381239b0595f98ccbb434b4305536dfdfd255c0fc96985

Comment

On 2025-03-18 11:42:16 Jacob Schneidewind wrote:

Notes:

- * Fixing link to H₂[17O] water preparation, specifying that water and Ag₂O are 10% 17O labelled
- * 9:25 deleting duplication
- * Putting dashes into empty fields (photos or observations etc.)
- * Including information on appearance of reactants (Ru(PNP) and Ag₂O etc.)
- * Making "light turned off" in grey
- * Specifying that stirring was turned off at 14:45
- * 14:55: fixing uM to um
- * 15:33: specifying room temperature water bath
- * at 14:55: writing work-up procedure
- * 13.03, 13:30: specifying that solid was transferred in air, and specifying that at this point the sample was weighed, and referring to AE-460-2 instead of "the solid"
- * Adding MeCN-d₃ to work-up reagents
- * Having heading for results
- * Adding future recommendation: using freshly prepared Ag₂O



Unique eLabID: 20250310-9e0bbe50a767eac9cc87147abd8cc955b206a975

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1857>