

Lauda - LOOP Thermoelectric thermostat, L100, 1

Date: 2025-09-29

Tags:

Lauda

Irradiation setup

EA

Category: Equipment

Created by: Ebrahim Abedini

Table

Device name	LOOP Thermoelectric thermostat
Model number	L100, S220011320
Paper, Information	Lauda operation manual.pdf
Location	CEEC II, lab E004 Irradiation set-up 1
Irradiation setup	Equipment - Irradiation setup, 1
Name in Instruments Calendar	/
Introduction	Alex, Nadja
Responsible Person	Alex, Nadja



Protocol on usage of the device	Protocol - Cleaning of double walled beaker
Additional Information	Fill it with 0.1 g/L Na₂CO₃ solution Remove and work gently with the inlet and outlet adapters. Connect the outlet of the Lauda to the lower connection of the double-walled beaker and the inlet to the upper one.
Status from	29/09/2025

Linked resources

Equipment - [Irradiation setup 1 \(CEEC II, E004\)](#)

Protocol - [Cleaning of double walled beaker](#)

Attached files

Lauda-L100-1.jpg

sha256: 366c0d540d3bdec5a2d1ce5159f9f1cd9aef7cc02a550413a080850779296403



Lauda operation manual.pdf

sha256: 9473256beea086f2b8c3ccd717390278dad446eb4bd58842398f699d33f5ceaa



Unique eLabID: 20250929-e29ef51ce5aa4a2236b499e39bd72d01882d4a0c

Link: <https://elab.water-splitting.org/database.php?mode=view&id=273>