

AE-361: 5 mM Stocksolution of KOH in water

Date: 2024-10-17

Tags: Degassed KOH Stocksolution AE

Category: Prep work

Status: Done

Created by: Alexander Eith

Reaction scheme/sample structure

Chemical Formula: H_2O
Molecular Weight: 18,01500

KOH

H_2O

Chemical Formula: HKO
Molecular Weight: 56,10530

Reagents

Name	CAS/Exp. Nr.	Amount [mmol]	Molar mass [g/mol]	Volume _{theo} [mL]	Volume [ml]	Density [g/mL]
KOH (1 M)	AE-243	0.1	56.1056	0.100	0.100	/
H2O (Mili-Q)	/	/	18.02	4.90	4.90	0.997

Procedure/observations

All steps were done according to [\[Protocol\] Schlenk Technique](#)

The flask was not dired prior to use

Date	Time	Step	Observations
17.10		A 10 mL Schlenk flask was set under argon.	
	11:35	To the flask water (not degassed) was added using 5 mL graduadet pipette	
		To the flask KOH (1 M) was added using a 100 μL Hamilton syringe	
	11:42 - 12:12	The solution was degassed according to Protocol - Degassing of fluids for 28 min	
		The flask was stored under argon in the fumehood	

Linked experiment

- [AE-243: 1 M Stocksolution of KOH in water](#)

Linked resources

Protocol - [Schlenk Technique](#)

Protocol - [Degassing of fluids in normal Schlenk flask](#)



Unique eLabID: 20241018-9e5ce94e58e39b41b519b0b659e1812a6ffe3481
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1430>