

AE-152: Degassing of MeCN

Date: 2023-09-11

Tags: Degassed AE

Status: Done

Created by: Alexander Eith

Reaction scheme/sample structure

Reagents

Name	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
Ethanol (SPS)	/	/	/	/	/	approx. 450 mL

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard Schlenk technique or in an argon filled glovebox (MBRAUN).

Date	Time	Step	Observations
11.09	15:50	A 500 mL-solvent flask was prepared.	
		Molar sieve 3 A (approx. 100 mL) was filled in the flask in argon counterstream	
		MeCN was collected from the SPS according to Protocol - Collecting of SPS solvents	
	16:55 - 17:35	Argon was sparged through the solution analog to [Protocol] Degassing of fluids for 40 min	

Analysis

Date	Time	Sample name	Analysis method	Solvent	Raw Data	Processed Data	Interpretation

Product characterization

Linked resources

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [KFT \(Karl-Fischer-Titration\)](#)

Protocol - [Collecting of SPS solvents](#)



Unique eLabID: 20230911-891565e76df8da057b7d9414460aa152dc2e5841
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=625>