

AE-426: Radiation of [Ru(CO)(OH)₂(PNP)] in water with liquid phase O₂ measurment

Date: 2025-01-23

Tags: Radiation O₂ PNP

[Ru(CO)(OH)₂(PNP)] AE liquidphase


Category: Two-photon

Status: Done



Created by: Alexander Eith

Reaction scheme/sample structure

Radiation of Dihydroxo



Chemical Formula: C₂₄H₄₅NO₃P₂Ru
Molecular Weight: 558,64552
[Ru(CO)(OH)₂(PNP^{tBu})]



Chemical Formula: C₂₄H₄₅NO₂P₂Ru
Molecular Weight: 542,64652
[Ru(CO)H(OH)(PNP^{tBu})]

Literature/reference experiments

Literature	https://doi.org/10.1039/D1EE01053K
Reproduction	/
Related experiments	Two-photon - AE-418: Radiation of [Ru(CO)(OH) ₂ (PNP)] in water with liquid phase O ₂ measurment Two-photon - AE-419: Radiation of water with liquid phase O ₂ measurment, analog to AE-418

Reagents

Name	CAS/Experiment number	Amount [μmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
[Ru(CO)(OH) ₂ (PNP{tBu})]	Two-photon - AE-382: Synthesis of [Ru(CO)(OH) ₂ (PNP ^{tBu})] using Ag ₂ O, 5 h (AE-382-3)	/	/	/	2.31	558.65	/
Milli-Q water	GD-018: Redegassing of H ₂ O	/	/	/	/	/	0.92

Irradiation Parameters

Power measured using [\[Power Meter\] 843-R-USB + 919P-020-12](#) unless specified otherwise.

	Light Source Name	Wavelength [nm]	Power Setting [mW]
First light source	Light Source - LCS-0365-76-22	365 nm	300
Second light source	Light Source - LCS-6500-65-22	white	200

Used beam combiner [Name or None]	Beam Combiner - LCS-BC25-0409
Irradiation distance [cm]	9
Thermostat temperature [°C]	25
Stirring speed [rpm]	1000
Irradiation start [s]	2880
Irradiation stop [s]	8720

Procedure/observations

All steps were carried out using [\[Protocol\] Schlenk Technique](#). All steps were carried out in the dark, unless other specified.

Date	Time	Step	Observations
23.01		A 10 mL screw cap schlenk flask was prepared	
	13:35	To the flask [Ru] was added under air	
		The flask was set under Ar	
	13:55	To the flask water (0.92 mL) was added	
		The mixture was stirred	brwonish suspension in water.jpg
	14:00	The flask stored in Equipment - Freezer Lab E004 CEEC II under Ar in the dark	

24.01		A screw cap 5 mL irradiation flask was prepared	
	9:30	The 10 mL flask was taken out of the freezer	
		The power of the LEDs 365 nm: 300 mW and white: 200 mW were adjusted according to [Protocol] Irradiation setup	power setup.jpg
	9:40	The 10 mL flask was stirred at 600 rpm	after stir.jpg no change
	10:00	The suspension was filtered into the 5 mL flask using a PA syringe filter (pore diameter = 0.2 µm)	yellowish solution after filt.jpg
	10:04	The Robust oxygen probe was attached to the flask using a bola fitting (diameter = 3 mm)	
	10:05	The flask was moved to the setup	in setup.jpg
	10:07	The logging was started (AE-426-1-Ch2)	
	10:08	The flask was flushed with argon for 130 s (strong Ar flow in gas phase)	
	10:11 - 10:35	The Schlenk was degassed according to Protocol - Degassing and O2 Measurement with a sensor cap with air-stable substances in the Photochem. setup	c(O2) = approx. -0.6 µmol/L during degas.jpg no change
	10:35	The degassing stopped by removing the PTFE cannula according to Protocol - Degassing and O2 Measurement with sensor cap with air stable substances	
	10:37	The Log (-1-Ch2) was stopped	c(O2) = approx. 0.2 µM AE-426-1.png AE-426-1-Ch2.txt
	10:37	The logging was started (AE-426-2-Ch2)	
	10:37 - 11:25	The sample was allowed to equilibrate according to Protocol - Degassing and O2 Measurement with sensor cap with air stable substances	Equilibration duration: 2880 s
	11:25	The light sources were turned on	
	13:03	The light source was turned off	Irradiation duration: 5840 s

	13:03 - 15:13	The sample was allowed to equilibrate according to Protocol - Degassing and O2 Measurement with sensor cap with air stable substances	Equilibration duration: 7834 s after irrad.jpg orange solution with small amount of solids
	13:54	The logging was stopped (-2-Ch2)	AE-426-2-Ch2.txt AE-426-2.png

Analysis

Date	Time	Sample name	Analysis method	Analytic device	Solvent	Raw Data	Processed Data	Comparative Data	Interpretation
24.01	10:07	AE-426-1-Ch2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	H2O	AE-426-1-Ch2.txt	AE-426-1.png	Two-photon - AE-402: Radiation of [Ru(CO)(OH)2(PNP)] in water with liquid phase O2 measurment (AE-402-1-Ch2 and AE-402-2-Ch2)	Degassing worked. After stop increase to 0.2 µM
24.01	10:37	AE-426-2-Ch2	Optical O2 detection	Equipment - Firesting Fiber-Optic Oxygen Meter 2 Channel	H2O	AE-426-2-Ch2.txt	AE-426-2.png	Two-photon - AE-402: Radiation of [Ru(CO)(OH)2(PNP)] in water with liquid phase O2 measurment (AE-402-3-Ch2) Two-photon - AE-418: Radiation of [Ru(CO)(OH)2(PNP)] in water with liquid phase O2 measurment (AE-418-2-Ch2)	Very similar to AE-402 and AE-418. Slow decrease during equilibration. Sharp drop to -0.6 at start of irrad, than graduale increase by 0.1 over approx. 5 min, than slow decrease similar to decrease before irrad. At stop of irrad drop to -0.05 µmol/L, than slow increase over time (maybe slowly flattening)

Result

Calibration using static degassed conditions did not help to get rid of the artefacts. Otherwise reproducible

Linked experiments

- [AE-236: Irradiation of PhPDA \(AE-234\), 1.5 mg/mg SDS, 2 mg/mL PhPDA](#)

- AE-237: Irradiation of PhPDA (AE-234), 1.5 mg/mg SDS, 2 mg/mL PhPDA I
- AE-245: Synthesis of poly(1,4-diphenylbutadiyne)
- AE-249: Irradiation of PhPDA (AE-245), 1.5 mg/mg SDS, 2 mg/mL PhPDA
- AE-254: Irradiation of PhPDA (AE-245), 1.5 mg/mg SDS, 2 mg/mL PhPDA reproduction
- AE-257: Synthesis of poly(1,4-diphenylbutadiyne)
- AE-258: Irradiation of PhPDA (AE-257), 1.5 mg/mg SDS, 2 mg/mL PhPDA
- AE-259: Irradiation of PhPDA (AE-257), 1.5 mg/mg SDS, 2 mg/mL PhPDA after 1 d
- AE-260: Irradiation of PhPDA (AE-257), 1.5 mg/mg SDS, 2 mg/mL PhPDA after 2 d
- AE-262: Irradiation of PhPDA (AE-257), 1.5 mg/mg SDS, 2 mg/mL PhPDA after 3 d
- AE-276: Synthesis of poly(1,4-diphenylbutadiyne), 1.2 g scale
- AE-279: Irradiation of PhPDA (AE-276-2), 1.5 mg/mg SDS, 2 mg/mL PhPDA
- AE-282: Irradiation of PhPDA (AE-276-1), 1.5 mg/mg SDS, 2 mg/mL PhPDA
- GD-018: Redegassing of H₂O
- AE-313: Synthesis of PPy via Fenton reaction with FeSO₄ and 10 eq. H₂O₂ I
- AE-332: Irradiation of PPy (AE-312), 2 mg/mL
- AE-333: Irradiation of PPy (AE-313), 2 mg/mL
- C3N4 - AE-398: Irradiation of C3N4-Co-0.2 from Lingli Ni
- C3N4 - AE-399: Irradiation of C3N4-Co-0.4 from Lingli Ni I
- Two-photon - AE-358: Radiation of [Ru(CO)(OH)₂(PNP{tBu})] (AE-310) 365 nm and white, 50 mM KOH
- Two-photon - AE-382: Synthesis of [Ru(CO)(OH)₂(PNPtBu)] using Ag₂O, 5 h
- Two-photon - AE-386: Radiation of [Ru(CO)(OH)₂(PNP{tBu})] (AE-382) 365 nm and white, 50 mM KOH in D₂O I
- Two-photon - AE-402: Radiation of [Ru(CO)(OH)₂(PNP)] in water with liquid phase O₂ measurment
- Two-photon - AE-418: Radiation of [Ru(CO)(OH)₂(PNP)] in water with liquid phase O₂ measurment

Two-photon - [AE-419: Radiation of water with liquid phase O2 measurment, analog to AE-418](#)

Linked resources

Archived - [Degassing procedure for photocatalytic water splitting by polymeric photocatalysts](#)

Beam Combiner - [Beam Combiner LCS-BC25-0409 -1](#)

Beam Combiner - [Beam Combiner LCS-BC25-0409 -2](#)

Equipment - [Firesting Fiber-Optic Oxygen Meter 2 Channel \(Firesting 1\)](#)

Equipment - [Freezer Lab E004 CEEC II](#)

Light Source - [UHP LED 405 nm](#)

Light Source - [UHP LED 6500 K \(white\)-1](#)

Light Source - [HP LED 405 -2](#)

Light Source - [HP LED 6500 K \(white\) -2](#)

Protocol - [Irradiation setup](#)

Protocol - [UV vis measurements \(BioTek PowerWave HT\)](#)

Protocol - [Determination of H2O2 by Spectroquant](#)

Protocol - [Degassing and O2 Measurement with sensor cap with air stable substances](#)

Protocol - [Determination of H2O2 in the suspensions of photocatalysts using H2O2 test strips](#)

Protocol - [Ultrasonication with the UP 200st from hilscher in ZAF 208](#)

Protocol - [H2O2 photometric detection with DPD-POD method for irradiated suspensions of photocatalysts](#)

Attached files

in-setup.jpg

sha256: 3320fdbcb65a3b9825626ad1e18a60e13bfa3a482938fd421ff19f01ab5532ee



after-irrad.jpg

sha256: 7bf0b8c784d4518ce3baf3c5b00fc103084bf05e3461120dbd96d9aeccbe18ea



after-stir.jpg

sha256: 9579fbd8275955944c6f060ef567d4b88d9543a3685c027ed7b5fac29374ea87



after-filt.jpg

sha256: 6c8103f1a03d189e044cb8c89ea4042a934ba51d8ac73cce2c4c439096202a42



during-degas.jpg

sha256: 5b79292cd11af7a9db74aa3528c2dc8287ff2d4d5425f12832693cd8127aca25



power-setup.jpg

sha256: ea1dff285ef1200f74af142eaf405ce3911865ad15bd5c99ef34680465f8941e




in-water.jpg

sha256: 20322ca5fa639168dff3812f4e67fe7e046e6620ef2bbe9d183c053eb2b780bc



AE-426-2.png

sha256: ed60260dc494a5ece03ee4822defda7c8c5f9e0c032a76c4e8ade60e2dc8cc6b




AE-426-2-Ch2.txt

sha256: 98f3de5ccd189002750cd0b2c7586acc06a2760a9a36ab26fbbeafc840c79adc

AE-426-1.png

sha256: e828f5a0529ee944f2c819d42094d21605175dacd2ea1f5a43f40bc345ff988c



AE-426-1-Ch2.txt

sha256: 89bc6165d86f9e535b9228595d4d88a1038c93cdf215b91362c9ce7a5dd2865d



Unique eLabID: 20250123-b835cba2958feb2c87aa2b7137e87e34b9643b12

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1664>