

GD-018: Redegassing of H2O

Date: 2024-06-11

Tags: Water GD Degassing

Status: Done

Created by: Gustav Döhring

Reaction scheme/sample structure

Literature/reference experiments

Literature	/
Reproduction	/
Related experiment	Experiment - GD-016: Redegassing of H2O

Reagents

Name	CAS Number / Experiment Number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]
H2O	Experiment - GD-016: Redegassing of H2O	/	/	/	/	/	/	aprox. 150

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard [Protocol - Schlenk Technique](#) or in an argon filled glovebox (MBRAUN).

Date	Time	Step	Observations
05.06	12:50-13:20	Argon was sparged through the solution analog to [Protocol] Degassing of fluids for 30 min	

Results

Degassing worked

Linked experiments

- [GD-003: Degassing of Milli-Q-Water](#)

- [GD-016: Redegassing of H2O](#)



Unique eLabID: 20240611-7c0fb21db64c654706f3ea6ad2d962a2f9a610fb
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=1094>