

KRA-104: Degassing of Milli-Q water

Date: 2024-01-15

Tags: Degassed solvent KRA

Status: Done

Created by: Kristína Rabatinová

Reagents

Name	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Volume [ml]
Milli-Q water	/	/	/	/	/	approx. 70 mL

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard Schlenk technique or in an argon filled glovebox (MBRAUN).

Date	Time	Step	Observations
15.01.24		A schlenk flask (100 mL) was prepared.	
		Milli-Q water was added (70 mL) to the schlenk flask.	
	10:50 - 11:30	Argon was sparged through the solution analog to [Protocol] Degassing of fluids for 20 min	

Linked experiment

- [AE-208: Degassing of Milli-Q water](#)

Linked resources

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [KFT \(Karl-Fischer-Titration\)](#)

Protocol - [Collecting of SPS solvents](#)



Unique eLabID: 20240115-bb92d8c796158de8135e77014575fa1f2dff3d63
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=802>