

# KRA-068: Degassing of DCM

Date: 2023-12-05

Tags: Degassed solvent KRA

Status: Done

Created by: Kristína Rabatinová

## Reagents

| Name | Amount [mmol] | Equivalents | Mass <sub>theo</sub> [mg] | Mass <sub>exp</sub> [mg] | Molar mass [g/mol] | Volume [ml] |
|------|---------------|-------------|---------------------------|--------------------------|--------------------|-------------|
| DCM  | /             | /           | /                         | /                        | /                  |             |

## Procedure/observations

All steps (unless stated otherwise) were carried out under argon using either standard Schlenk technique or in an argon filled glovebox (MBRAUN).

| Date     | Time          | Step                                                                                       | Observations |
|----------|---------------|--------------------------------------------------------------------------------------------|--------------|
| 05.12.23 | 15:30         | J.Young schlenk flask was prepared                                                         |              |
|          | 16:00         | Molar sieve was added to the flask (approx. 1/4)                                           |              |
| 06.12.23 | 14:00         | Collecting at SPS according to <a href="#">Protocol - Collecting of SPS solvents</a>       |              |
|          | 15:05 - 15:40 | The DCM was degassed according to <a href="#">Protocol - Degassing of organic solvents</a> |              |

## Linked resources

Protocol - [Degassing of organic solvents in Young-type flask](#)

Protocol - [KFT \(Karl-Fischer-Titration\)](#)

Protocol - [Collecting of SPS solvents](#)



Unique eLabID: 20231206-4039c3d3684071bd1250c2550d80bebfba7209a8  
Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=747>