

AE-617: Radiation of [Ru(CO)(OH)2(PNP{tBu})] (AE-610), kinetic

Date: 2025-11-11

Tags: Dihydroxo Radiation PNP
P(tBu)N(py)P(tBu) [Ru(CO)(OH)2(PNP)]
NMR AE 1H 31P Irradiation
photocatalysis

Category: Two-photon

Status: Done


Created by: Alexander Eith

Objective

kinetic study of formation of hydride and side products during irradi of [Ru(CO)(OH)2(PNP{tBu})]

Reaction scheme/sample structure

Radiation of Dihydroxo



ChemDraw File (linked): [AE-590.cdxml](#)

Literature/reference experiments

Literature	https://doi.org/10.1039/D1EE01053K
Reproduction	/
Similar experiments	Two-photon - AE-543: Radiation of [Ru(CO)(OH)2(PNP{tBu})] (AE-456-1), kinetic, max. conc Two-photon - AE-597: Radiation of [Ru(CO)(OH)2(PNP{tBu})] (AE-456-1), kinetic, max. conc

Reagents

Name	CAS Number / Experiment Number	Inventory number	Amount [mmol]	Equivalents	Mass _{theo} [mg]	Mass _{exp} [mg]	Molar mass [g/mol]	Density (g/ml)	Volume [ml]	Concentration [mM]	Obtained concentration [mM]
[Ru(CO)(OH)2(PNP{tBu})] (max. conc approx. 1.44 mM) in water	Two-photon - AE-610: Stocksolution of [Ru(CO)(OH)2(PNP{tBu})] in water	/	/	/	/	/	558.6455	/	0.60	max. conc (approx. 1.44 mM)	approx. 1.44

Irradiation Parameters

Power measurement was performed in experiment [Two-photon - AE-616: Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\] \(AE-610\), kinetic](#)

	Name
Used Set-up	Equipment - Advanced irradiation setup V1.0 I
Irradiation set-up number	Equipment - Irradiation setup 2 (CEEC II, E004)

	Light Source Name	Used power source	Wavelength [nm]	Power Setting [mW]	Analogue Setting [-]
First light source	Light Source - UHP LED 365 nm-2	Power Sources - BLS-18000-1 1	365	500	9.42
Second light source	Light Source - UHP LED 470 nm-2	Power Sources - BLS-13000-1E	470	500	6.97

Used beam combiner [Name or None]	Beam Combiner - Beam Combiner LCS-BC25-0409 -1
Irradiation distance [cm]	6.5
Thermostat temperature [°C]	25
Stirring speed [rpm]	/
Irradiation duration	1, 3, 20, 24 h

Procedure/observations

All steps (unless stated otherwise) were carried out under argon using standard [Protocol - Schlenk Technique](#)

The synthesis/reactions were performed in the dark by either wrapping the reaction flasks in aluminum foil or turning the lab light and light in the fume hood off, when the flasks were uncovered.

The NMR tube was placed in the back holder of the holder for 2 NMR tubes

Date	Time	Step	Observations	Pictures
-------------	-------------	-------------	---------------------	-----------------

06.11	14:55	A Young type NMR was prepared according to Protocol - Preparation of NMR Sample (from step 3) using 0.60 mL of AE-610 , a Prep work - AE-489: Capillary production containing D2O and H3PO4	AE-617	/
	15:20	The NMR sample an NMR was measured	AE-617-0H slightly yellow solution	AE-617-0H.jpg
10.11	09:05	The sample was placed in the setup using the Photoreactor - AE-558: Design and print of holder and holder disk for 2 NMR tubes Holder_2NMR_tubes	no change	before irradiation 1.jpg
	09:38	The irradiation was started	/	during irradiation 1.jpg
	10:38	The irradiation was stopped after 1h	slightly more orange	after irradiation 1.jpg
	09:20	An NMR was measured	AE-617-1H	AE-617-1H.jpg
11.11	09:44	The sample was placed in the setup using the Photoreactor - AE-558: Design and print of holder and holder disk for 2 NMR tubes Holder_2NMR_tubes	no change	before irradiation 2.jpg
	09:44	The irradiation was started	/	during irradiation 2.jpg
	12:44	The irradiation was stopped after 3h	more orange	after irradiation 2.jpg
	12:58	An NMR was measured	AE-617-4H	AE-617-4H.jpg
12.11	12:51	The sample was placed in the setup using the Photoreactor - AE-558: Design and print of holder and holder disk for 2 NMR tubes Holder_2NMR_tubes	no change	before irradiation 3.jpg
	12:51	The irradiation was started	/	during irradiation 3.jpg
13.11	08:51	The irradiation was stopped after 20h	no big change	after irradiation 3.jpg
	09:01	An NMR was measured	AE-617-24H	AE-617-24H.jpg
	17:17	The sample was placed in the setup using the Photoreactor - AE-558: Design and print of holder and holder disk for 2 NMR tubes Holder_2NMR_tubes	no change	/
	17:18	The irradiation was started	/	/

14.11	17:21	The irradiation was stopped after 24h	no big change, maybe slightly darker	/
	17:39	An NMR was measured	AE-617-48H , no big change, maybe slightly darker	AE-617-48H.jpg

Analysis

Date	Time	Sample name	Analysis method	Analytical device	Solvent	Raw Data	Python script	Processed Data	Comparative Data	Interpretation
06.11	21:12	AE-617-0H	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H2O	AE-617-0H.zip	/	AE-617-0H_10.nmrium	/	As expected, before irradiation, dihydroxo is main species
10.11	14:22	AE-617-1H	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H2O	AE-617-1H.zip	/	AE-617-0H_10.nmrium	/	No big difference
11.11	13:32	AE-617-4H	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H2O	AE-617-4H.zip	/	AE-617-0H_10.nmrium	/	Hydride observed
13.11	12:32	AE-617-24H	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H2O	AE-617-24H.zip	/	AE-617-0H_10.nmrium	/	Significant hydride formation
17.11	11:37	AE-617-48H	¹ H, ³¹ P{ ¹ H}quant NMR	IAAC 400 MHz	H2O	AE-617-48H.zip	/	AE-617-48H_10.nmrium	/	Full conversion

Sample	Yield (hydride species, cis) [%]	Yield (hydride species, trans) [%]	Yield (dihydroxo) [%]	Yield (main side product, 57 ppm) [%]	Yield (soluble side products) [%]	Yield (other side products) [%]	Data
AE-617-1H	0	0	93	0	4	3	AE-617_Yield_Determination.xlsx
AE-617-4H	2	0	62	3	7	27	AE-617_Yield_Determination.xlsx
AE-617-24H	17	0	11	32	12	27	AE-617_Yield_Determination.xlsx
AE-617-48H	22	0	0	33	7	38	AE-617_Yield_Determination.xlsx

Results

Worked, significant product after 4 h, full conversion after 48 h, similar yields compared to [Two-photon - AE-616: Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\] \(AE-610\), kinetic](#).

2 differences: higher conversion here after 4 h, otherwise conversion analog to AE-616

In AE-616 decrease in main side product from 24 to 48 h, here nearly constant

Linked experiments

Prep work - [AE-489: Capillary production](#)

Prep work - [AE-549: Degassing of Milli-Q water](#)

Two-photon - [AE-597: Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\] \(AE-570-1\), kinetic, max. conc](#)

Two-photon - [AE-603: Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\] \(AE-570-1\), kinetic, max. conc, reproduction](#)

Two-photon - [AE-610: Stocksolution of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\] in water](#)

Two-photon - [AE-616: Radiation of \[Ru\(CO\)\(OH\)2\(PNP{tBu}\)\] \(AE-610\), kinetic](#)

Linked resources

Beam Combiner - [Beam Combiner LCS-BC25-0409 -1](#)

Equipment - [Advanced irradiation setup V1.0 I](#)

Equipment - [Advanced power measurment setup V1.0 I](#)

Equipment - [Irradiation setup 2 \(CEEC II, E004\)](#)

Light Source - [UHP LED 365 nm-2](#)

Light Source - [UHP LED 470 nm-2](#)

Power Meter - [843-R-USB + 919P-020-12](#)

Power Sources - [BLS-13000-1E](#)

Power Sources - [BLS-18000-1 1](#)

Protocol - [Schlenk Technique](#)

Attached files

AE-617-48H.jpg

sha256: 5608d0a6e059a3437b5971e457acda5a30ac7d3f2683690a8a98aaf07d6d555b



AE-617_Yield_Determination.xlsx

sha256: ee8fd8c65a4fb53d87f2b706504282d015f57278114e4e648633eca47601bf8f

AE-617-48H_10.nmrium

sha256: 39a4d22c8cec0fecf90eba55d4d5c9a6cd55a25f71f8a5d5662b411ecac4c08f

AE-617-48H.zip

sha256: f1d5fa516d6ddfaccce57e6afbfc1dc0524749c340c98169103e3bb87307cd51

AE-617_Yield_Determination.xlsx

sha256: 3a748bc886311e2c8b764f3ce37d61d405609479f57228c090a4276d70d573b8

AE-617-0H_10.nmrium

sha256: 15a6d2b7e2fd28111ce2c724590348bdb88abb8b6eaf3df23845b689bad6c297

AE-617-24H.zip

sha256: 11762b6c2c725a82e5acd7708d20879aee54f0396807585eba930bf17cb66eb9

AE-617-4H.zip

sha256: 190aca21bd3ff5155c6cce75d2204e9abf1ce650b35c42133a782c7ee0a43722

AE-617-24H.jpg

sha256: c668cb9415c2175681cdde80eaf50072547c099c8080e6f38a124f89fd1b7bb9



after irrad 3.jpg

sha256: 796ed31e21f1c652ee5151f3885eb4cf80fe2d5470b8636fa8a7d151468a0de1



during irrad 3.jpg

sha256: 86ae82664ce8e0649397d13ff34f02eef63e469c2897aee2a2ac26aff9b1c933



before irradi 3.jpg

sha256: ac867c24203dd684fbddbe81e5120e73d41ca8d058b086cf29eecea3875dd966



AE-617-1H.zip


sha256: 57c5f6b74b759a00c4bb0456b610216552c2f34a5534ba140bf4fa9bdd179534

AE-617-0H.zip

sha256: c2eeda4747cabbcc00c64d0ab9954ae3c9e13ff879a336d5f413bf410a631484


AE-617-4H.jpg

sha256: 2e54046e488175fb0dbe8d1ba7dc6e3939baf15102eddd5ef6c59b7ce35779e



after irradi 2.jpg

sha256: bcf9d9b4b79d7cd7246454b51d108a80f350683e93e662ffbbbbb60ef14c1a4b3



during irradi 2.jpg

sha256: f909ad3d917421ce65b7bbf30777db4fc14a70d905dfdc042e9c4344d3489791



before irradi 2.jpg

sha256: 39798cb2bf4725106e5d2d79cca9573ce3a0f0333f231d97e7d0ebb6dd2cc653



AE-617-1H.jpg

sha256: b1b69ec1873ce18a9dbe21de4a1649349fbe7088f87a610d999413feacfed87d



after irrad 1.jpg

sha256: 87ee2e52e3e6a348f3e1d5dd2182e0aac65242f4b573e04e6c8b7c915f415832




during irrad 1.jpg

sha256: bcf3e754aca8de66a0545e41ebc133a72d3bf9b7d5f4a3b05e2e8000c0670961




before irrad 1.jpg

sha256: d1521b5d2f6d0b470d7076c3e42ea5643f03b43b3498d6e95b4e02792b35bcd9



AE-617-0H.jpg

sha256: 08a23002fa7d5d10f22445b3c929cd8de699f21da5155678ed3fabe95d4eeade



AE-590(9).cdxml

sha256: 2bb17a23d6ac9133244c0395480ff7b2d05f24ff8c1a41893e7170ba2481ccdc



Unique eLabID: 20251111-23c061d46b47458f1787bd225557cae381749cde

Link: <https://elab.water-splitting.org/experiments.php?mode=view&id=3413>